BioMonitor 2

ProMRI®

Cardiac Monitor | with Fast Insert Tools | BIOTRONIK Home Monitoring

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Table of Contents

Product Description	2
Intended Medical Use	2
System Overview	3
Diagnostic Functions	5
General Safety Instructions	6
Operating Conditions	6
Possible Complications	7
Possible Risks	7
Handling	8
Insertion Procedure	8
Precautionary Measures while Programming	12
Follow-Up	12
Patient Information	12
Replacement Indications	13
Explantation and Device Replacement	13
Disposal	13
Parameters	14
Arrhythmia Detection Parameters	14
Home Monitoring Parameter Settings	15
Technical Data	16
Mechanical Characteristics	16
Electrical Characteristics	16
Battery Data	17
Legend for the Label	18

1 Product Description

Intended Medical Use

Intended Use

BioMonitor 2-AF and BioMonitor 2-S are the names of implantable cardiac monitors for the monitoring and automatic recording of the heart rhythm. Their primary purpose is to provide early detection and diagnostics of symptoms of arrhythmias, such as atrial tachycardias and the causes of syncope, which can be clinically manifested.

Note: BioMonitor 2 does not have a pacing function.

The FIT accessory kit consists of two tools: Fast Insert Tool 1 and 2 (FIT 1 and FIT 2). They are used for forming the device pocket and for subsequent positioning of the cardiac monitor in the subcutaneous left pectoral area. The use of these tools ensures an optimal anatomical implantation site, which is a requirement for recording meaningful subcutaneous ECGs.

The accessory is an integral component of the medical product and may be used exclusively for the insertion of BioMonitor 2-AF and 2-S. It can be ordered separately as a sterile back-up kit.

Form of diagnosis

The heart rhythm is continuously and automatically recorded and monitored. Possible detection types are atrial tachycardias (only BioMonitor 2-AF), high ventricular rate, asystole, bradycardia, and sudden ventricular rate drop. Depending on the preset parameters, subcutaneous ECGs and other data may be recorded.

The patient can also trigger the recording of subcutaneous ECGs using the

Remote Assistant accessory if subjectively symptomatic episodes occur.

The recordings are transmitted to the BIOTRONIK Home Monitoring Service Center.

This enables physicians to perform diagnosis management at any time, day or night.

Required expertise

In addition to having basic medical knowledge, the user must be thoroughly familiar with the operation of the implantable cardiac monitor and the specified conditions for its use. Only qualified medical specialists having the special knowledge are permitted to implant the BioMonitor 2 and make diagnoses.

Guidelines of cardiological societies

It is recommended that the indications published by the DGK (German Cardiac Society) and the ESC (European Society of Cardiology) are observed. This also applies to the guidelines published by the HRS (Heart Rhythm Society), the ACC (American College of Cardiology), the AHA (American Heart Association), and other national cardiology associations.

Indications

Generally approved differential diagnostics methods, indications, and recommendations for the medical use of implantable cardiac monitors apply to BioMonitor 2-AF and 2-S including the FIT 1 and FIT 2 accessories.

The BioMonitor 2-AF and 2-S cardiac monitors are implantable monitoring systems that record subcutaneous ECGs. Recording is activated both automatically and by the patient. Their use is indicated for the following cases:

- Clinical symptoms or increased risk of cardiac arrhythmias
- Temporary, unexplained symptoms that may indicate cardiac arrhythmia

Contraindications

There are no known contraindications.

However, the particular patient's state of health determines whether a subcutaneous device will be tolerated long-term.

System Overview

Device family

BioMonitor 2-AF (atrial fibrillation) and BioMonitor 2-S (syncope) are implantable cardiac monitors, which differ in their diagnostic range of functions: BioMonitor 2-AF can also sense and record episodes of atrial tachycardias.

These devices are not available in every country.

Parts

The system consists of the following parts:

- Device with flexible lead body
- Accessories FIT 1 and FIT 2
- Programmer and current software version for the device
- Optionally the Remote Assistant accessory for triggering recordings by the patient


Illustration of tools and the device

- FIT 1 (left): Tool (handle and insertion area) for forming the device pocket
- FIT 2 (middle): Tool for controlled insertion and positioning of the device
- BioMonitor 2 (right): Device


Note: The tools in the two differently packaged accessory sets vary, as can be seen in the following two figures. Please also see the possibility of distinguishing the accessory sets using the sterile order number in section "Sterile packaging of accessories FIT 1 and FIT 2".

- Please also refer to the quick reference guide included in the package contents.

Tool (individually packaged in sterile bags) and device:

Tool (packaged together in a blister) and device. This accessory set can be optically identified by a thickened, ergonomically designed handle of the FIT 1.

Note: Please also see the different handling when inserting the BioMonitor 2: "Place device on FIT 2 and press if necessary".

Cardiac monitor

The device consists of a solid housing part and a flexible part, the lead body. This flexible lead body is composed of silicone and carries the antenna for Home Monitoring in addition to the lead.

The housing of the device has the approximate shape and size of a USB flash memory stick and is made of biocompatible titanium, welded from the outside and thus hermetically sealed. It is coated in silicone. At the rounded end of the housing there is an opening in the coating, so that the metal housing forms the antipole to the lead tip. The device has an overall length of 8.8 cm, which is approximately congruent with the sensing vector and correlates linearly with the sensing amplitude.

Programmer

Implantation and follow-up are performed with a portable BIOTRONIK programmer using the current PSW software version 1503.A or higher. During the implantation, the standard program is activated in the device on initial programming via the programmer. The programmer is used to set parameter combinations, as well as for interrogation and saving of data from the device. Subcutaneous ECG, markers, and functions are displayed simultaneously on the color display.

Note: The programmer's ECG display must not be used for diagnostics, because it does not meet all the standard requirements for diagnostic ECG devices (IEC 60601-2-25).

Telemetry

Telemetric communication between the device and the programmer can be carried out by applying a programming head (PGH).

BIOTRONIK Home Monitoring®

The BIOTRONIK cardiac monitors provide a complete diagnosis management system:

- With Home Monitoring, diagnostic information as well as device technical data are automatically and wirelessly sent to a stationary or mobile transmitter via an antenna in the lead body. The data are encrypted and sent from the transmitter to the BIOTRONIK Home Monitoring Service Center (HMSC) through a cellular phone network.
- The received data are deciphered and evaluated. Each physician can set the criteria for evaluation to be used for each patient and can configure the time of notification via e-mail, SMS or fax.
- A clear overview of the results of this analysis is displayed for the attending physicians on the protected Internet platform Home Monitoring Service Center (HMSC).
- Data transmission from the device is performed with a daily device message.
- Device messages that indicate special events in the heart or in the device are forwarded at the preset time.
- A test message can be initiated at any time using the programmer to immediately check the Home Monitoring function.

BioMonitor 2 order numbers

Not all device types are available in all countries.

Silicone housing coating:

Device	Order number
BioMonitor 2-AF	398493
BioMonitor 2-S	398494

FIT 1 + 2 order number

Accessory kit as a reserve

Accessories	Order number
Fast Insert Tools FIT 1 + FIT 2	408906

Note: Please also see the possibility of distinguishing the accessory sets using "Illustration of tools and the device" and "Sterile packaging of accessories FIT 1 and FIT 2".

Package contents

The storage package includes the following:

- Sterile packaging with device
- Sterile packaging with accessory (FIT 1 and FIT 2 either individually in a sterile bag or both together in a single sterile blister)
- Serial number label
- Patient ID card
- Quick reference guide for the insertion of BioMonitor 2 using the accessories

Note: The technical manual pertaining to the cardiac monitor and accessories is either included in hard copy form in the storage package or is available in digital form on the Internet.

Each sterile packaging includes the following:

- Device
- FIT 1 accessory
- FIT 2 accessory

Diagnostic Functions

General overview

- Automatic functions make it easy and fast to insert, configure, and check the BioMonitor 2-AF and 2-S cardiac monitors.

Detection and data storage

- The parameters for the sensing settings, such as high pass filter, target sensing threshold, or noise window, can be adjusted individually for each patient.
- The sensing parameters are combined in one program (SensingConsult). In addition to the possible individual program, there are preconfigured programs for PVC detection, for T-wave suppression, and for management of variable amplitudes and a setting named Standard.
- The signals are automatically recorded and stored once a detection type is activated and detection occurs.
- Multiple detection types can be activated simultaneously.
- A total of 55 individual episodes with a length of at least 40 s each can be stored automatically. The maximum storage period for an individual episode is 60 s. A total of 4 recordings triggered by the patient with a duration of at least 7.5 min can be stored. The recording includes 7 min of pre-episode history and 0.5 min of post-episode history relative to the time of triggering. The maximum recording duration for an individual episode is 10 min.
- The device can store episodes with subcutaneous ECGs with a minimum overall length of 60 min.
- When performing in-office follow-ups using the programmer, the subcutaneous ECG is indicated with markers after applying the programming head.

Range of functions

The two device types differ in terms of their range of functions:

Diagnostic functions	BioMonitor 2-AF	BioMonitor 2-S
Longest AF episode	x	—
AF burden		
AF trend		
AF time distribution		
AF duration		
Ventricular rate during AF		
Atrial fibrillation	x	—
High ventricular rate	x	x
Bradycardia	x	x
Sudden rate drop	x	x
Asystole	x	x
Patient trigger	x	x

Home Monitoring functions

Important medical information includes the following:

- Sustained atrial arrhythmias (only BioMonitor 2-AF)
 - Sustained ventricular arrhythmias
 - Current statistics
 - Periodically recorded subcutaneous ECGs that are transmitted according to an individually adjustable timing interval in addition to the regular device message
- This is a necessary condition for performing Home Monitoring-supported follow-ups.

2 General Safety Instructions

Operating Conditions

Technical manuals

The following technical manuals provide information about usage of the device systems:

- Technical manual for the device
- Technical manual for the HMSC
- Technical manuals for the programmer
- Technical manuals for the user interface
- Technical manuals for cables, adapters and accessories
- Technical manuals are either included in hard copy form in the product package or are available in digital form on the Internet: manuals.biotronik.com
- Follow all relevant technical manuals.
- Keep technical manuals for later use.

Note: Observe the included quick reference guide.

Care during shipping and storage

- Devices must not be stored or transported close to magnets or sources of electromagnetic interference.
- Take the effects of storage period and service time for the device into consideration, see Battery data.
- Observe the use by date of the accessory kit.

Temperature

Extremely low and high temperatures affect the service time of the device's battery. The temperature data also apply to the FIT 1 and FIT 2 accessories.

- Permitted for shipping and storage:
-10°C to +45°C (-40°F to 158°F)

Sterile delivery

The device and the accessories have been gas-sterilized (ethylene oxide). Sterility is guaranteed if the blister/sterile bags and quality control seals have not been damaged.

Device sterile packaging

The device is double packaged in 2 separately sealed sterile blisters. The inner blister is sterile on the outside so that it can be transferred and remain sterile during implantation.

Sterile packaging of accessories FIT 1 and FIT 2

The accessories are available in two versions:

- FIT 1 and FIT 2 are both individually mounted to a transport insert and sealed in their own sterile bag (sterile order numbers 408999 and 409000).
- Or FIT 1 and FIT 2 are packaged together in a single sterile blister (sterile order number 429000).

Single use only

The device and the FIT 1 and FIT 2 accessories are intended for single use only.

- Do not use the device and accessories when one of the packages is damaged.
- The device and accessories must not be resterilized and reused.

Possible Complications

General information on medical complications

Complications for patients and device systems generally recognized among practitioners also apply to BIOTRONIK devices.

- Known complications are foreign body rejection phenomena, local tissue reactions, migration of the device, accumulation of fluid in the device pocket, transdermal erosion, and/or infections. Primary sources of complication information include current scientific and technological knowledge.
- It is impossible to guarantee the efficacy of implantable cardiac monitors.

Possible technical failures

Technical failure of a device system cannot be entirely ruled out. Possible causes can include the following:

- Component failure of the device or an accessory
- Battery depletion

Electromagnetic interference (EMI)

Any device can be sensitive to interference, for example, when external signals are sensed as intrinsic rhythm.

- BIOTRONIK devices have been designed so that their susceptibility to EMI is minimal.
- Under unfavorable conditions, especially during diagnostic and therapeutic procedures, sources of interference may induce such a high level of energy into the device that the data recording can be influenced or the device may be damaged.

Possible Risks

Procedures to avoid

The following procedures must be avoided as they may cause harm to the patient or damage the device and, as a result, put the system functionality at risk:

- Therapeutic ultrasound
- Transcutaneous electrical nerve stimulation
- Hyperbaric oxygen therapy
- Applied pressures higher than normal pressure

Risky therapeutic and diagnostic procedures

If energy from an external source is introduced to the body for diagnostic or therapeutic purposes, then the device can be subjected to interference and its function can be impaired.

Harmful effects may occur, for example, during electrocautery, HF ablation, HF surgery, shortwave and microwave therapy, or lithotripsy. Influences on the device cannot always be detected immediately.

If risky procedures cannot be avoided, the following should be observed at all times:

- Electrically insulate the patient.
- Do not approximate energy nearby the device.
- Monitor the patient during and after every intervention.

External defibrillation

The device is protected against the energy that is normally induced by external defibrillation. Nevertheless, any implanted device may be damaged by external defibrillation. The sensing properties may change as a result.

- Position the adhesive electrode anterior-posterior or vertical to the connection line from the device to the heart as well as minimum 10 cm away from the device.

Radiation therapy

The use of radiation therapy must be avoided due to possible damage to the device and the resulting impaired functional safety. If this type of therapy is to be used anyway, prior risk/benefit analysis is absolutely necessary. The complexity of influencing factors such as different sources of radiation, a variety of devices and therapy conditions makes it impossible to issue directives that guarantee radiation therapy without an impact on the device. The EN 45502 standard pertaining to active implantable medical devices requires the following information in combination with therapeutic ionizing radiation:

- Adhere to instructions for risky therapeutic and diagnostic procedures.
- Shield device against radiation.
- After applying radiation, double-check the device system to make sure it is functioning properly.

Note: Please contact BIOTRONIK with questions during the risk/benefit analysis.

Magnetic resonance imaging

Magnetic resonance imaging must be avoided due to the associated high-frequency fields and magnetic flux density. This could result in damage or destruction of the device by strong magnetic interaction and damage to the patient by excessive warming of the body tissue in the area surrounding the device.

Under certain conditions and while maintaining mandatory measures, magnetic resonance imaging can be performed to protect the patient and device.

- The ProMRI® manual – MR conditional device systems – contains detailed information on safely conducting an MR scan.
 - Download the digital manual from the Internet: manuals.biotronik.com
 - Order the printed manual at BIOTRONIK.
- Does approval as "MR conditional" apply in your country or region?
Ask for current information at BIOTRONIK.

3 Handling

Insertion Procedure

Having parts ready

The following parts that correspond to the requirements of EC Directive 90/385/EEC are required:

- BIOTRONIK device
- BIOTRONIK tools for forming the device pocket and for subsequent insertion and positioning of the device
- BIOTRONIK programmer and approved cables
- External multi-channel ECG device
- Keep spare parts for all sterile components.

Unpacking the device

WARNING

Inadequate function due to defective device

If a device still in packaging is dropped on a hard surface during handling, electronic parts could be damaged.

- Use a replacement device.
 - Return the damaged device to BIOTRONIK.
-
- Peel the sealing paper off of the outer blister at the marked position in the direction of the arrow. The inner blister must not come into contact with persons who have not sterilized their hands or gloves, nor with non-sterile instruments!
 - Take hold of the inner blister by the gripping tab and take it out of the outer blister.
 - Peel the sealing paper off of the sterile inner blister at the marked position in the direction of the arrow.

Unpacking accessories

Unpack tools from each of the 2 sterile bags: FIT 1 and FIT 2 are each mounted on a transport insert that has been sealed in a sterile bag.


- Open the sterile bag at the marked position. The transport insert and the tools must not come into contact with persons who have not sterilized their hands or gloves or with non-sterile instruments!
- Have the transport insert removed by a person who has sterilized his/her hands or gloves.
- Remove FIT 1 and FIT 2 from the transport insert.

Unpack tools from the blister: FIT 1 and FIT 2 have been individually packaged in a sterile blister.

- Open blister at the marked position. The tools must not come into contact with persons who have not sterilized their hands or gloves or with non-sterile instruments!
- Remove FIT 1 and FIT 2 from the blister.

Implantation site

The cardiac monitor is inserted subcutaneously on the left side of the patient. There are 2 primary positions depending on the patient's anatomy:

Position **A**: device pocket at an angle of about 45° to the sternum, parallel to the cardiac vector

Position **B**: device pocket in the left parasternal region


Note: In exceptional cases, the position beneath the left chest can be selected.

Sequence of the procedure

Follow these steps to insert the cardiac monitor; please also refer to the included quick reference guide:


Make an incision for the device pocket

- Using a scalpel, cut an incision at least 15 mm wide for the device pocket.


Hold FIT 1 correctly

- Hold the tool by the handle, between your thumb and index finger.


Form the device pocket

- Insert the tool subcutaneously up to the handle to form the device pocket.


- Do not form any skin folds with the other hand, so that the device pocket does not deform upon removal.


Place device on FIT 2 and press if necessary

- Please also observe the quick reference guide to distinguish handling!
- Handling of FIT 2 from the sterile bag (sterile order numbers 408999 and 409000): Place the device onto FIT 2 in the preformed opening.
- Handling of FIT 2 from the blister (sterile order number 429000): Place the device into the preformed recess on the FIT 2 and clip it in by pressing the flexible lead body in the guide beneath the undercut.


Position FIT 2 with the device

- Hold the device and FIT 2 securely between your thumb and index finger.


Insert FIT 2 with the device

- Insert the device on the tool into the device pocket, at least up to the circumferential welded seam on the device; the flexible lead body must lie on the tool during insertion.


- While pushing the implant into the device pocket, make sure it does not slip over the tool; the flexible lead body should not bend and must rest on the tool for the duration of the insertion.

Insert the device completely


- Put your thumb on the still protruding end of the device and completely insert it.

Withdraw FIT 2 from the device pocket

- Press slightly on the skin above the end of the device to fixate it. Then withdraw FIT 2 from the device pocket.

Make sure that the device remains in the device pocket and does not slip.

End the procedure

- Check the signals.
- Close the device pocket.

Applying the programming head

Put programming head (PGH) on and turn to align if necessary to ensure correct telemetry.

- Make sure the PGH is positioned correctly.

Establishing telemetry contact

- The device requires some time for interrogation once the programming head is applied.

Activate diagnostics

- Confirm the standard program is loaded on the programmer.

Precautionary Measures while Programming

Note when adjusting bigeminy rejection

For sensitive AF detection, **bigeminy rejection** is preset to the parameter **Standard**. If bigeminy rejection is set to the parameter **Aggressive**, the sensitivity may be slightly reduced depending on the patient.

- Adjust the setting for this parameter based on the current patient data.

Follow-Up

Follow-up intervals

Follow-ups must be performed at regular, agreed upon intervals.

- Following the ingrowth phase, approximately 3 months after insertion, the first follow-up should be carried out by the physician using the programmer (in-office follow-up).
- The next in-office follow-up should be carried out once a year and no later than 12 months after the last in-office follow-up.

Follow-up with BIOTRONIK Home Monitoring®

Monitoring using the Home Monitoring function does not serve to replace regular in-office appointments with the physician required for other medical reasons.

Home Monitoring-supported follow-up can be used to functionally replace in-office follow-up under the following conditions:

- The patient was informed that the physician must be contacted if symptoms worsen or if new symptoms arise despite the use of the Home Monitoring function.
- Device messages are transmitted regularly.
- The physician decides whether the data transmitted via Home Monitoring with regard to the patient's clinical condition as well as the technical state of the device are sufficient. If not, an in-office follow-up has to be carried out.

Possible early detection due to information gained via Home Monitoring may necessitate an additional in-office follow-up. For example, the data may provide early indications of a foreseeable end of service (EOS). Furthermore, the data could provide indications of previously unrecognized arrhythmias.

Follow-up with the programmer

Use the following procedure for in-office follow-up:

1	Record and evaluate the ECG.
2	Interrogate the device.
3	Check the sensing function.
4	Evaluate status and, if necessary, statistics and subcutaneous ECGs.
5	If necessary, customize the program functions and parameters and transmit the edited program to the device.
6	Print and document follow-up data (print report).
7	Finish the follow-up for this patient.

Patient Information

Patient ID card

A patient ID card is included in package contents.

- Provide the patient with the patient ID card.
- Request that patients contact the physician in case of uncertainties.

Prohibitive signs

Premises with prohibitive signs must be avoided.

- Educate the patient on prohibitive signs.

Possible sources of interference

Electromagnetic interference should be avoided in daily activities. Sources of interference should not be brought into close proximity with the device.

- Educate the patient on sources of electromagnetic interference which include special household appliances, safety locks, anti-theft alarm systems, cell phones and transmitters.
- Request patients to do the following:
 - Use cell phones on the side of the body that is opposite of the device.
 - Keep the cell phone at least 15 cm away from the device both during use and during storage.

Recording via Remote Assistant

Patients can also trigger recordings manually using the Remote Assistant accessory. They need to know the following:

- The physician must explain to patients how to handle the Remote Assistant based on the technical manual.
- When can the Remote Assistant be used?
In case of symptoms such as severe dizziness, indisposition, or palpitations, as well as after fainting.
- How is a recording triggered?
The technical manual for the Remote Assistant describes the handling and the interpretation of the LED signals.
- What happens when a recording is triggered?
The subcutaneous ECG of the last 7.5 min is stored. This includes 7 min of pre-episode history and 0.5 min of post-episode history relative to the time of triggering. Between 40 s and 60 s are transmitted at the daily transmission time via Home Monitoring.
- What has to happen after triggering a recording?
The patient has to contact the physician for a follow-up.

Replacement Indications

Possible battery levels

- BOS: Beginning of Service
- ERI: Elective Replacement Indication
- EOS: End of Service

Elective Replacement Indication (ERI)

Home Monitoring is still available for 2 weeks once the device enters ERI. After this period, the device stops transmitting messages to the Home Monitoring Service Center. The recording of statistics and episodes is also stopped after entering ERI.

EOS replacement indication

The end of service has been reached.

Explantation and Device Replacement

Explantation

- Interrogate the device status.
- Remove the device using state-of-the-art technology.
- Explanted devices are biologically contaminated and must be disposed of safely due to risk of infection.

Device replacement

Basic principles:

- The device must not be resterilized and reused.

Cremation

Devices should not be cremated.

- Explant the device before the cremation of a deceased patient.

Disposal

Disposal of the device

BIOTRONIK takes back used products for the purpose of environmentally safe disposal.

- Clean the explant with a solution of at least 1% sodium hypochlorite.
- Rinse off with water.
- Send the cleaned explant to BIOTRONIK.

Disposal of the tools

Accessories FIT 1 and FIT 2:

- Used accessories must be disposed of as medical waste in an environmentally friendly and proper manner.

4 Parameters

Arrhythmia Detection Parameters

Atrial fibrillation detection settings

Valid for BioMonitor 2-AF, the following can be set:

Parameter	Range of values	Standard	Factory
Atrial fibrillation (AF)	ON; OFF	ON	OFF
AF sensitivity	Low; Medium; High; In case of change of the AF expert parameters: individual	Medium	Medium

Additional AF expert parameters can be set individually:

Parameter	Range of values	Standard	Factory
RR variability limit	6.25; 12.5; 18.75	12.5% (medium)	
Onset/resolution window	8/16; 16/24; 24/32	8/16 cycles (medium)	
Onset intervals	5 ... [2] ... 23	5 cycles (medium)	
Resolution intervals	1 ... [2] ... 7	1 cycle (medium)	
Confirmation time	1 ... [1] ... 6; 10; 20; 30	6 min (medium)	
Bigeminy rejection	OFF; Standard; Aggressive	Standard (medium)	

The AF expert parameters are preset as follows:

Parameter	Low	Medium	High
RR variability limit	18.75%	12.5%	6.25%
Onset/resolution window	16/24 cycles	8/16 cycles	8/16 cycles
Onset intervals	9 cycles	5 cycles	5 cycles
Resolution intervals	3 cycles	1 cycle	1 cycle
Confirmation time	6 min	6 min	6 min
Bigeminy rejection	Standard	Standard	Standard

High ventricular rate

The following can be set:

Parameter	Range of values	Standard	Factory
High ventricular rate (HVR)	ON; OFF	ON	OFF
HVR limit	110 ... [10] ... 200	180 bpm	180 bpm
HVR counter	8 ... [4] ... 24; 32	16 cycles	16 cycles

Bradycardia

The following can be set:

Parameter	Range of values	Standard	Factory
Brady zone limit	OFF; 30 ... [5] ... 80	40 bpm	OFF
Brady duration	5 ... [5] ... 30	10 s	10 s

Sudden rate drop

The following can be set:

Parameter	Range of values	Standard	Factory
Sudden rate drop (SRD)	OFF; Rate decrease: 20 ... [10] ... 70%	OFF	OFF
SRD sensitivity	Low; Medium; High; In case of change of the SRD expert parameters: individual	Medium	Medium

Additional SRD expert parameters can be set:

Parameter	Range of values	Standard	Factory
Baseline intervals	48; 64; 128; 256	64 cycles (medium)	
Rate-drop intervals	8; 16; 32	16 cycles (medium)	

The SRD expert parameters are preset as follows:

Parameter	Low	Medium	High
Baseline intervals	256 cycles	64 cycles	48 cycles
Rate-drop intervals	32 cycles	16 cycles	8 cycles

Asystole duration

The following can be set:

Parameter	Range of values	Standard	Factory
Asystole duration	OFF; 2 ... (1) ... 10	3 s	OFF

Patient trigger

The following can be set:

Parameter	Range of values	Standard	Factory
Patient trigger	ON; OFF	ON	OFF

Resting rate period

The following settings apply to all detection types:

Parameter	Range of values	Standard	Factory
Start resting period	00:00 ... (1:00) ... 23:00	2:00 hh:mm	
Duration of resting period	0.5 ... (0.5) ... 12	4.0 h	

Sensing settings: SensingConsult

Individual programs (SensingConsult) with expert sensing parameters can be set as follows:

Parameter	Range of values	Standard	Factory
Recordings & real-time	Full morphology; Sensing signal	Full morphology	
Interference interval	100 ... (10) ... 200	100 ms	

Parameter	Range of values	Standard	Factory
Sensing settings (Sensing-Consult)	Standard; Variable amplitude; PVC detection; T-wave suppression; In case of change of the expert sensing parameters: individual	PVC detection	
Target sensing threshold	25; 35; 40; 50	25% (PVC detection)	
Sensing high-pass filter	10; 18; 24	10 Hz (PVC detection)	

The expert sensing parameters are preset as follows:

Parameter	Programs			
	Standard	Variable amplitude	PVC detection	T-wave suppression
Target sensing threshold	50%	35%	25%	40%
Sensing high-pass filter	10 Hz	10 Hz	10 Hz	18 Hz
Interference interval	100 ms	100 ms	100 ms	100 ms
Recordings & real-time	Full morphology	Full morphology	Full morphology	Full morphology

Home Monitoring Parameter Settings

General settings

The following settings apply to all detection types:

Parameter	Range of values	Standard	Factory
Home Monitoring	ON; OFF	ON	OFF
Time of transmission	STD 00:00 ... (00:30) ... 23:30 hh:mm	STD	STD
Periodic subcutaneous ECG for Home Monitoring-supported follow-up	OFF; Cycle duration: Selection; 1; 7; 30; 60; 90; 120; 180 Selection: 1st date ... 5th date	30 days	OFF
		30 days	OFF

HM episode trigger

Episode triggers for Home Monitoring can be set for all detection types:

Parameter	Range of values	Standard	Factory
Atrial fibrillation (AF)	ON; OFF	ON	OFF
High ventricular rate (HVR)	ON; OFF	ON	OFF
Bradycardia	ON; OFF	ON	OFF
Sudden rate drop (SRD)	ON; OFF	OFF	OFF
Asystole	ON; OFF	ON	OFF
Patient trigger	ON; OFF	ON	OFF

Note: When detection type AF, HVR or Bradycardia is set, the oldest, the most recent, and the longest episodes are saved.

When detection type Asystole or Sudden rate drop is set, the oldest and the two most recent episodes are saved.

When Patient trigger is set, the oldest and the three most recent episodes are saved.

5 Technical Data

Mechanical Characteristics

Housing dimensions

Device	L x H x D [mm]	Volume [cm ³]	Mass [g]
BioMonitor 2-AF/2-S	88.4 x 15.2 x 6.2	5	10.1

Accessory dimensions

FIT 1 and FIT 2 packed in separate sterile bags:

Product	L x H x D [mm]
FIT 1	179 x 15 x 3
FIT 2	128.2 x 17 x 1.5

FIT 1 and FIT 2 packed in one blister:

Product	L x H x D [mm]
FIT 1	203 x 15 x 3
FIT 2	129.8 x 17 x 2.9

X-ray identification

A0

Materials in contact with body tissue

- Housing: titanium
- Lead body: silicone
- Lead tip: titanium
- Housing coating: silicone
- Additional components (adhesive): silicone
- FIT 1 and FIT 2: plastic (PA/PE)

Electrical Characteristics

Components and input values

Electrical characteristics determined at 37°C, 500 Ω:

Circuit	Hybrid electronics with VLSI-CMOS chip
Input impedance	> 50 kΩ

Housing shape

The housing has an elongated, flattened shape similar to a USB flash memory stick, connected to a narrow, flexible lead body.

Electrically conductive surfaces

There is an electrically conductive surface on the housing as well as on the lead body:

Uncoated area of the electrically conductive housing surface [mm ²]	Uncoated area of the electrically conductive lead tip [mm ²]
53	53

Telemetry

Telemetry data for Home Monitoring:

MICS frequencies	Maximum power of transmission
402 – 405 MHz	< 25 µW -16 dBm

International radio certification

Devices with BIOTRONIK Home Monitoring® are equipped with an antenna for wireless communication.

- Telemetry information for Australia:

This product is in compliance with the Australian "Radiocommunications Act 1992" and therefore it is labeled according to the "Radiocommunications (Compliance Labeling - Devices) Notice".

- Telemetry information for Japan:

In accordance with Japanese law, this device has been assigned an identification number under the "Ordinance concerning certification of conformity with technical regulations etc. of specified radio equipment", Article 2-1-8.

R: 202-SMC062

- Telemetry information for the USA:

This transmitter is authorized under the Medical Device Radiocommunication Service (in part 95 of the FCC Rules) and must not cause harmful interference to stations operating in the 400.150 to 406.000 MHz band in the Meteorological Aids (i.e., transmitters and receivers used to communicate weather data), the Meteorological Satellite, or the Earth Exploration Satellite Services and must accept interference that may be caused by such stations, including interference that may cause undesired operation. This transmitter shall be used only in accordance with the FCC Rules governing the Medical Device Radiocommunication Service. Analog and digital voice communications are prohibited. Although this transmitter has been approved by the Federal Communications Commission, there is no guarantee that it will not receive interference or that any particular transmission from this transmitter will be free from interference.

This device will be registered with Federal Communications Commission under the following number:

FCC ID: QR1BM2

Battery Data

Battery characteristics

The following data is provided by the manufacturers:

Manufacturer	LITRONIK GmbH 01796 Pirna, Germany
Battery type	LiS 2460
System	LiMnO ₂
Device type	BioMonitor 2-AF and 2-S
Battery voltage at BOS	3.1 V
Open-circuit voltage	3.2 V
Nominal capacity	440 mAh
Remaining capacity at ERI	40 mAh
Usable capacity until EOS	400 mAh

Power consumption

- Power consumption with factory settings: 3.2 µA
- Power consumption in transmission mode: 3 µA
- Power consumption from implantation to ERI: 6.1 µA
- Power consumption after ERI: 4.65 µA

Average service time

Based on the technical specifications of the battery manufacturer and the device-specific detection parameters, the average service time of BioMonitor 2-AF and 2-S is 4 years until reaching ERI. This calculation is based on sensing at 100 bpm and daily transmission of one device message via Home Monitoring.

The remaining life of BioMonitor 2-AF and 2-S is 2 months after entering ERI and before reaching EOS. Home Monitoring is still available for 2 weeks after this time. ERI is displayed on the programmer during follow-up and transmitted via Home Monitoring. All statistics and recording functions are deactivated.

Shortening of the service time after long storage period

Depending on the storage period, the service time from the beginning of service BOS to the replacement time ERI decreases as follows:

- After 1 year: by 4 months
- After 1.5 years: by 6 months

Legend for the Label

Meaning of the symbols

The label icons symbolize the following:

	Manufacturing date		Use by
	Storage temperature	REF	Order number
SN	Serial number	PID	Product identification number
CE	CE mark	LOT	Lot number
	Contents		Follow the instructions for use!
manuals.biotronik.com 		Follow the electronically available instructions for use!	

	Manufacturer	USA Distributor: Distributor in the USA
Rx only	Caution: Federal (U.S.A.) law restricts this device to sale by, or on the order of, a physician.	

STERILE	Sterilized with ethylene oxide
----------------	--------------------------------

	Do not resterilize		Single use only. Do not re-use!
	Do not use if packaging is damaged		Non-sterile

	Transmitter with non-ionizing radiation at designated frequency
	MR conditional: Patients can be examined using an MR scan under defined conditions if the device system in their body consists of devices whose packaging is labeled with this symbol.

TP2	Compatibility with telemetry protocol version 2 of BIOTRONIK Home Monitoring
------------	--

	Implantable cardiac monitor		
	FIT 1: Accessory for preparing the device pocket		FIT 2: Accessory for inserting the device
	Single sterile blister		

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All rights reserved. Specifications subject
- ® All product names in use may be trademarks or
registered trademarks held by BIOTRONIK