



Product Overview

Dental Implant Systems | Biomaterials



SICace®



SICmax®



**SICmax
onepiece**



SIC invent Around the World

What began more than 10 years ago with a small group of doctors and entrepreneurs is today a company group with a total of 26 subsidiaries and distribution partners worldwide.

Implants from Implantologists.





SIC invent



Dear Partners and Friends,

More than 10 years ago, my father and I founded SIC invent. What began with a Swiss corporation became a globally active company. The concept for success is based on our company slogan "Implants from Implantologists".

The SIC "Schilli Implantology Circle" was founded at the same time as the SIC invent group. SIC is an internationally organised network of opinion leaders and users of the system, jointly responsible for a powerful SIC invent product portfolio and its reliable application on patients, based on a worldwide continuous training concept.

Virtually all system components, which you will find in this product overview, have been developed in broad agreement with the SIC "Schilli Implantology Circle". Before they are presented in this catalogue, they have obtained the necessary evidence via studies and non-interventional studies at leading universities and in the practices and clinics of SIC members.

At this point, I would like to express my sincere gratitude to all SIC members. Without this "Think Tank" and user-oriented advisory group, a premium product and training portfolio would not be possible.

Thanks to our members and fellows amongst others, we will launch MY-SIC online this autumn. MY-SIC is a web-based platform which not only intends to promote continuing education of dental surgeons, prosthodontists and dental technicians but also as communication centre for professional exchange with colleagues.

You have in front of you our newly designed company catalogue in which we present not only the familiar range of products but also numerous innovations.

The individual approach

Let us inspire and convince you, now and in the future, with our products innovations and optimisations.

At this stage I would like to thank you for the confidence you have placed in us, and I look forward to continuing our cooperation.

Sincerely
Georg Schilli

A handwritten signature in black ink, appearing to read 'G. Schilli'.

President of the Supervisory Board
Chief Executive Officer
SIC invent AG

Schilli Implantology Circle Cooperation

Dear Partners and Friends,

Implantology has changed dentistry. It is, like the surgical joint replacement, a product of modern bone surgery. In 1958, this was completely restructured by the Arbeitsgemeinschaft für Osteosynthese (AO = Association for the Study of Osteosynthesis). Their rules also apply to us: the biomechanical principles of function provide orientation, the vitality of the bone must be preserved. The aim is the most atraumatic surgical procedure possible.

Incorporation of an implant is a biological process, and we must provide the requirements for it to take place smoothly. Every detail of our approach is therefore important. The procedure is optimised and errors are avoided by providing clear surgical protocols. But despite all the schematic optimisation every case remains an individual case. This in particular applies to the subsequent prosthetic restoration. When prosthodontics cooperates with surgery, prosthodontics dominates as it determines the function and aesthetics.


The implant system plays a major role in this process which is influenced by the individual factors of the patient. It is a standard product not only from a technical point of view but also for legal reasons. In this case, continuous optimisation is also a prerequisite for lasting success.

This is one of the tasks of the SIC "Schilli Implantology Circle". The international, interdisciplinary discussion forum integrates current theoretical and scientific research results in the implementation of practice-relevant systems and synchronises them with practical experience.

Another function of the SIC "Schilli Implantology Circle" is also to pass on this success in the form of continuous training.

We are, as are the members of the circle, very enthusiastic about implantology and want to improve the process and product through discussions with competent specialists. Disadvantages and faults as well as possible improvements and prospects are discussed at regional and an international level and scientific investigations and checks are initiated. Organisation and trust are therefore prerequisites for ensuring that everyone profits from this approach and that we can provide our patients with even more reliable and better help.

Sincerely
Prof. Dr. Wilfried Schilli



Founding Member Schilli Implantology Circle



SIC invent Overview Implant Systems	p. 8
SIC invent Surgical Guidelines	p. 12
SIC invent Prosthetic Guideline	p. 16
SIC invent Implants	p. 18
SICace®	p. 20
SICmax®	p. 21
SIC Cover Screws	p. 22
SICmax onepiece	p. 22
SIC Provisional Implants	p. 24
SIC invent Prosthetics	p. 26
SIC Gingiva Shapers	p. 28
SIC Gingiva Shapers Ø 3.3 mm – cylindrical	p. 28
SIC Gingiva Shapers Ø 3.3 mm – concave, anterior	p. 28
SIC Gingiva Shapers Ø 3.3 mm – concave, posterior	p. 28
SIC Gingiva Shapers Ø 4.2 mm – cylindrical	p. 29
SIC Gingiva Shapers Ø 4.2 mm – concave, anterior	p. 29
SIC Gingiva Shapers Ø 4.2 mm – concave, posterior	p. 29
Impression Technique	p. 30
SIC Impression Technique Ø 3.3 mm, reposition	p. 30
SIC Impression Technique Ø 4.2 mm, reposition	p. 30
Accessories – Transfer Caps	p. 30
SIC Impression Technique Ø 3.3 mm, Open Tray Technique	p. 31
SIC Impression Technique Ø 4.2 mm, Open Tray Technique	p. 31
Fixation Screws for Open Tray Technique	p. 32
SIC Insertion Tool	p. 32
SIC Lab Implants	p. 32
Abutments	p. 33
SIC P2F Abutment Dr. Galip Gurel	p. 33
SIC Standard Abutments Ø 3.3 mm, anterior, straight	p. 34
SIC Standard Abutments Ø 3.3 mm, anterior, 15° angle	p. 34
SIC Standard Abutments Ø 3.3 mm, anterior, 25° angle	p. 35
SIC Standard Abutments Ø 3.3 mm, posterior, straight	p. 35
SIC Standard Abutments Ø 3.3 mm, posterior, 15° angle	p. 36
SIC Standard Abutments Ø 4.2 mm, anterior, straight	p. 36
SIC Standard Abutments Ø 4.2 mm, anterior, 15° angle	p. 37
SIC Standard Abutments Ø 4.2 mm, anterior, 25° angle	p. 37
SIC Standard Abutments Ø 4.2 mm, posterior, straight	p. 38
SIC Standard Abutments Ø 4.2 mm, posterior, 15° angle	p. 38
SIC “Flex Star” Universal Cast-to Abutment	p. 39
SIC “Flex Star” Universal Cast-to Abutment	p. 39
SIC “White Star” Zirconium Oxide Abutment	p. 40
Accessories for Abutments	p. 40
SIC Multi-Unit Abutment System “Safe on Four”	p. 41
SIC Bar and Bridge Abutment “Safe on Four”	p. 42
SIC Standard Abutment “Safe on Four”	p. 42
SIC Accessories for Multi-Unit Abutment System	p. 43
SIC Crown Base “Safe on Four”	p. 43
SIC Bar and Bridge Abutment	p. 45
SIC Crown Bases and Bar Coping	p. 46
SIC Fixation Screws	p. 46
SIC Soldering Implant	p. 46
SIC Bars	p. 46

CAD/CAM	p. 47
SIC Bonding Bases	p. 47
SIC Wax-Up Bases	p. 48
SIC Scan Adapters	p. 49
Accessories for CAD/CAM	p. 49
Locator®	p. 50
Locator® Attachments	p. 50
Locator® Replacement Male	p. 51
Locator® Accessories	p. 51
Ball and Socket Attachments	p. 52
Inner Matrices	p. 52
Accessories Ball and Socket Attachments	p. 53
Steco Magnetic Attachments	p. 53
SICmax onepiece Accessories	p. 54
SIC Accessories Provisional Implants	p. 55
SIC Prosthetic Tools	p. 56
SIC invent Standard Surgery	p. 58
Surgical Tray	p. 60
Drill Systems	p. 60
SIC Marking Drill	p. 60
SIC Short Drills	p. 61
SIC Drills for Depth Stop	p. 61
SIC Depth Stops for Drills up to Ø 3.25 mm	p. 62
SIC Depth Stops for Drills up to Ø 4.25 mm	p. 62
SIC Brackets for Depth Stops	p. 63
SIC Crestal Drills	p. 63
SIC Reamers	p. 64
SIC Bone Taps	p. 64
SIC Bone Condenser	p. 65
SIC Surgical Accessories	p. 66
SIC Insertion Tools Surgery	p. 68
SIC invent Guided Surgery	p. 70
SIC Surgical Tray, Guided Surgery	p. 72
GS Gingiva Punches	p. 72
Drill Systems	p. 73
GS Countersinks	p. 73
GS Pilot Drills	p. 73
GS Extension Drills	p. 74
GS Bone Taps	p. 75
GS Insertion Tools	p. 75
GS Drill Keys	p. 76
GS Sleeves	p. 76
SIC invent Biomaterials	p. 78
Bone Grafting	
SIC nature graft	p. 79
Certificates	p. 80

SIC invent Implant Systems

SICace®

■ SICace®

The SICace® threaded cylindrical implant is characterised by an innovative, self-tapping screw design which follows the applied and recommended values of the "Arbeitsgemeinschaft für Osteosynthese" (AO - Association for the Study of Osteosynthesis).

It is suitable for all indications in oral implantology, even difficult anatomical situations, and for use in conjunction with augmentation techniques. The SIC drill system guarantees atraumatic preparation of the implant site and, together with the thread morphology, it enables high primary stability of the implant to be achieved, regardless of bone quality. Use of the SICace® implant can be recommended unreservedly in bone of D1 to D3 quality. The microstructure and degree of purity of the "SICmatrix" SIC surface ensure secure and lasting osseointegration.

The basic cylindrical shape of the implant has a conical taper in the crestal region which generates integrated "platform switching".

The prosthesis is attached via an internal precision hexagon. The small degrees of freedom of this abutment connection in combination with long parallel-walled guide surfaces, ensure maximum stability of the implant-abutment interface and provide lasting protection against screw loosening.




This is the all-round implant with outstanding clinical long-term results

Microstructure of the SIC surface – "SICmatrix" – in combination with Grade 4 titanium guarantees reliable, long-lasting osseointegration

Self-tapping screw design for all indications in oral implantology

- SIC drill system for an atraumatic preparation of the implant site
 - Basic cylindrical shape with apical conical taper for easy insertion of the implant
- Integrated "platform switching" for convenient prosthetic handling**
- Inner precision hex with long guide surfaces for maximum stability of the implant-abutment interface and a screw connection protected against continuous loading
 - Flexible and precise prosthetic components for all indications





■ SICmax®

When the SICmax® threaded cylindrical implant design was developed, priority was given to its suitability for “soft bone”. It is therefore to be used preferably for bone qualities D2 to D4. The implant site is prepared with the conventional SIC drill system.

A further feature is the greatly rounded implant tip without a direct thread cut. The SICmax® implant system is therefore particularly suitable for use in the upper posterior region especially with all forms of sinus lift.

The basic shape of the implant is cylindrical as regards to the thread flanks and the upper part was designed with conical “thickening” of the threaded core with a micro thread. In conjunction with the thread morphology, this guarantees very great primary stability. As a result, the implant is also suitable for immediate implantation.

A decompression zone was incorporated in the crestal emergence area which passes over into the integrated “platform switching”. These are construction details that improve the preservation and long-term stability of the peri-implant tissue. Like the SICace® implant, the implant body has the proven “SICmatrix” surface and an identical high-precision internal hexagon connection. The SICmax® implant can be restored prosthetically with all components of the SICace® system.

The solution in soft bone

- Implant specially designed for use in “soft bone”. It is therefore to be used preferably for bone qualities D2 to D4
- The implant site is prepared just as the SICace® implant (identical surgery tray, see SIC one4all Concept™)
- **Basic cylindrical shape of the implant with crestal micro thread guarantees very great primary stability. As a result, the implant is also suitable for immediate implantation**
- **Greatly rounded implant tip without a direct thread cut for use in the upper posterior region, especially with all forms of sinus lift**
- Integrated “platform switching” for convenient prosthetic handling
- Internal precision hexagon with long parallel-walled guide surfaces for maximum stability of the implant-abutment interface and a screw connection protected against continuous loading
- Microstructure of the SIC surface – “SICmatrix” – in combination with Grade 4 titanium guarantees reliable, long-lasting osseointegration
- Restorations are feasible with all prosthetic components of the SICace® implant system, flexible and precise solutions for all indications



SIC invent Implant Systems

SICmax onepiece

■ SICmax onepiece

The SICmax onepiece threaded cylindrical implant is a reduced diameter addition to the SIC implant systems.


The SICmax onepiece implant system allows economical and minimally invasive restorations by means of an O-ring attachment or the cementable SIC Retention Attachment. In addition, the multifunctional prosthetic interface in conjunction with the range of other prosthetic components enables restoration with fixed cemented bridges or single tooth restorations when the available space is reduced. The implant can be used especially with bone qualities D1 to D3, and the SIC "SICmatrix" surface guarantees a lasting stable implant-bone bond.

The basic form of the implant is cylindrical, and it has a "thickening" of the threaded core with a micro thread in the crestal region. The tulip-shaped gingival emergence passes over into a conical functional part which ends as a sphere to accept O-ring attachments.

The precise, anti-rotation prosthetic connection allows secure cementing of all available prosthetic superstructures and is optimised for the demands of CAD/CAM technology.




One-piece implant – Unique prosthetic flexibility

- 
- Multifunctional interface, ball connection and retention area provided by torx connection and conical fit for cemented abutments
 - Transgingival region with a height of 2 mm
 - Reinforced thread core with microthread for cortical penetration
 - Cylindrical main thread with large surface and aggressive apex
 - Microstructure of the SIC surface – "SICmatrix" – in combination with Grade 4 titanium guarantees reliable, long-lasting osseointegration

■ SIC Provisional Implants

The concept for temporary immediate prosthetic restoration is to avoid undesirable loading of simultaneously implanted definitive implants and to increase patient comfort.

Patients' expectations of implant-based restorations are extremely high nowadays. Temporary edentulousness is not accepted and removable temporary restorations are often refused. The SIC provisional implants provide a treatment concept that allows an economical and fixed interim restoration. The SIC Provisional Implants are available in a diameter of 2.3 mm and lengths of 10 mm and 13 mm. Both an O-ring attachment for fixing a prosthesis and a tapered abutment for cementing a bridge are available for prosthetic restoration. The only instruments required are a drill and an insertion instrument.



SIC invent Surgical Guidelines

SICace®

SICmax®

	<div> <div>3.4</div> <div>3.7</div> </div>	<div> <div>4.0</div> <div>4.2</div> </div>	<div> <div>4.5</div> <div>4.7</div> </div>	<div> <div>5.0</div> <div>5.2</div> </div>
SICace® SICmax®				
Standard Drills	6.0 to 14.5	6.0 to 14.5	6.0 to 14.5	6.0 to 14.5
Short Drills	6.0 to 11.5	6.0 to 11.5	6.0 to 11.5	6.0 to 11.5
Pilot Drill 2.0	x	x	x	x
Ext. Drill 2.8	x	x	x	x
Ext. Drill 3.1	x	x	x	x
Crestal Drill 3.3	x			
Bone Tap 3.4	(x)			
Ext. Drill 3.25		x	x	x
Crestal Drill 3.75		x		
Bone Tap 4.0		(x)		
Ext. Drill 3.75			x	x
Crestal Drill 4.25			x	
Bone Tap 4.5			(x)	
Ext. Drill 4.25				x
Crestal Drill 4.75				x
Bone Tap 5.0				(x)

(x) = optional, according to bone quality

	<div> <div>4.0</div> </div>	<div> <div>4.5</div> </div>	<div> <div>5.0</div> </div>
SICace® 6.0 mm Short Implants			
Standard Drills	6.0 to 14.5	6.0 to 14.5	6.0 to 14.5
Short Drills	6.0 to 11.5	6.0 to 11.5	6.0 to 11.5
Pilot Drill 2.0	x	x	x
Ext. Drill 2.8	x	x	x
Ext. Drill 3.1	x		
Ext. Drill 3.25		x	x
Reamer 4.0	x		
Ext. Drill 3.75			x
Reamer 4.5		x	
Reamer 5.0			x

Short and long drills with depth stops





1 Marking Drilling:
After detaching the mucosa, the implant site is spot drilled using the marking drill. The drilling depth is variable and is maximum the length of the implant.
Speed max. 800 rpm.



2 Pilot Drilling:
Using the Ø 2.0 mm pilot drill, the implant length and axial alignment are then determined. The drilling depth can be checked optically using the depth markings on the drills or using the optional depth stop system. Speed max. 800 rpm.




3 Depth Measurement:
Check of the preparation depth and axial alignment using the Ø 2.0 mm depth gauge.




4 Initial Extension Drilling:
The pilot drill hole is extended using the "Smart Drill" Ø 2.8 mm extension drill. The axial alignment can still be slightly adjusted at this stage. Speed max. 800 rpm.



5 Axial Alignment:
Check of the axial alignment using the Ø 2.8 mm depth gauge.



6 Further Extension Drilling:
The use of the respective SIC Ø 3.1 / 3.25 / 4.25 mm extension drills depends on the planned implant diameter (see overview top left). Speed max. 800 rpm. (We recommend a lower speed with each of the last extension drills. Below 60 rpm external cooling using a chilled, sterile, physiological saline solution (NaCl) or ringer solution is no longer required.)



7 Crestal Drill:
The depth of drilling depends on the bone quality (see graphic). Speed max. 650 rpm. (We recommend a lower speed when using the crestal drill. Below 60 rpm external cooling using a chilled, sterile, physiological saline solution (NaCl) or ringer solution is no longer required.)





8 Bone Tap:
An initial thread should be cut in cases with highly cortical bone structures (D1/D2). The depth of the thread cut is normally half the length of the implant. The thread can be tapped manually or using a handpiece. Speed max. 35 rpm.



9 Implant Placement
The implant insertion is performed with the insertion tool for angle piece, direct (max. 25 rpm) or additionally with the TR insertion tool and torque ratchet until it is flush with the bone level.



10 Implant closure:
Finally the implant is sealed using the cover screw and saliva-proof wound closure.





Align buccally

Caution:
The final position of the implant should always be aligned so that a corner of the outer hex of the insertion tool and implant insertion post faces buccally.


SIC invent Surgical Guidelines

SIC Provisional Implants



1 Pilot Drilling:

After making a flap and marking the implant location, the Ø 1.5 mm pilot drill is to be used to establish the implant length and alignment. Especially for dense bone (D1), a 1.7 mm pilot drill is useful. Max. speed 800 rpm.



2 Provisional Implant Insertion:

Inserting the implant can be done machine driven with the Insertion Tool for Angle Piece directly (max. 25 rpm) or additionally with the SIC TR Adapter for Angle Piece for Ratchet. The implant should be inserted until the crestal end of the thread.



3 Prosthetic Restoration:

After insertion, the provisional implant can be combined either with an O-Ring Attachment or a Conical Abutment. (Both are useful in case of overdenture prosthesis and bridgework).



1 Pilot Drilling:

After making a flap and marking the implant location, the Ø 2.0 mm pilot drill is to be used to establish the implant length and alignment.




2 Extension Drill optional:

Optional: With highly cortical bone structures (D1/D2), the first drilling step should also be drilled using Ø 2.5 mm of the extension drill "Smart Drill" Ø 2.8 mm. In case of very dense bone, the final diameter of 2.8 mm can also be used to a maximum of the first depth groove of the drill at 7.5 mm. Max. speed 800 rpm.



3 Insertion:

The implant insertion is performed with the insertion tool for angle piece, direct (max. 25 rpm) or additionally with the TR insertion tool and torque ratchet until the bottom edge of the instrument has reached the desired gingival level. The optimum gingival niveau and insertion depth are to be defined preoperative and dependent on the planned prosthetic solution.



4 Prosthetic Restoration:

There are a variety of options for prosthetic restoration available from O-Ring Attachment for fixation of fixed-removable restorations to cementable abutments or even customized CAD/CAM superstructures.

Align
buccally



Caution:
The final position of the implant should always be aligned so that the clip of the insertion tool faces buccally.

SIC invent Prosthetic Guideline



Gingiva Shaper

- cylindrical
- concave anterior
- concave posterior



Impression Technique

Reposition/Open Tray Technique




SIC P2F
Abutment
Dr. Galip Gurel
p. 33

**Standard
Abutments**
p. 34


**CAD/CAM
Abutments**
p. 47

**"Flex Star"
Abutments**
p. 39

**“Flex Star”
(bridge)
Abutments**
p. 39



Single tooth restoration



Bridgework



**Complete Resoration
non-removable**



**Complete Restoration
removable**

Tightening Torque	
Gingiva shaper	5 Ncm (hand tight)
Impression post	5 Ncm (hand tight)
All Fixation Screws and Attachments	20 Ncm
SIC Fixation Post "Safe on Four" with hex 2.4 mm	30 Ncm

**Multi-Unit
Abutment
System
"Safe on Four"**
p. 41




**Bar and Bridge
Abutments**
p. 45



**"White Star"
Abutments**
p. 40





**Locator[®]
Attachments**
p. 50



**Steco
Magnetic
Attachment**
p. 53



**Ball and
Socket
Attachments**
p. 52



Implants



SICace®	p. 20
SICmax®	p. 21
SIC Cover Screws	p. 22
SICmax onepiece	p. 22
SICmax onepiece Accessories	p. 23
SIC Provisional Implants	p. 24
SIC Provisional Implants Accessories	p. 25









SIC one4all Concept™

All three SIC invent implant lines are compatible with the instruments of the SIC Surgical Tray.
With this "SIC one4all Concept™", we meet the demands of the dental practitioner for ease of use and cost-effectiveness.

Implants

SICace®






SICace Screw Implants incl. Cover Screw

Diameter	 Ø 3.4 mm	 Ø 4.0 mm	 Ø 4.5 mm	 Ø 5.0 mm
Material: Titanium Grade 4				
Ø Prosthetic connection [mm]	3.3		4.2	
Geometry of the connection	Inner Hex 2.3			
Length [mm]	Art. No.			
6.0*		935140 [NEW]	935141 [NEW]	935142 [NEW]
7.5	935182	935183	935150	935184
9.5	935170	935174	935151	935178
11.5	935171	935175	935152	935179
13.0	935172	935176	935153	935180
14.5	935173	935177	935154	935181







*Please consider the differing drilling protocol on page 12.

SICmax®




SICmax Screw Implants incl. Cover Screw

Diameter	 Ø 3.7 mm	 Ø 4.2 mm	 Ø 4.7 mm	 Ø 5.2 mm
Material: Titanium Grade 4				
Ø Prosthetic connection [mm]	3.3		4.2	
Geometry of the connection	Inner Hex 2.3			
Length [mm]	Art. No.			
6.0		935265 [NEW]	935266 [NEW]	935267 [NEW]
7.5	935270	935275	935285	935280
9.5	935271	935276	935286	935281
11.5	935272	935277	935287	935282
13.0	935273	935278	935288	935283
14.5	935274	935279	935289	935284





Cover Screws

SIC “Augmentation Plus” Cover Screw		incl. SIC “Augmentation Plus” Fixation Screw		
Compatible with	 SICace Ø 3.4 mm SICmax Ø 3.7 mm	 SICace Ø 4.0 mm SICmax Ø 4.2 mm	 SICace Ø 4.5 mm SICmax Ø 4.7 mm	 SICace Ø 5.0 mm SICmax Ø 5.2 mm
Material: Titanium Grade 5				
Ø Implant connection [mm]	3.3		4.2	
Diameter [mm]	3.3		4.2	
Total height [mm]	2.4			
Geometry of the connection [mm]	cylindrical 2.25			
	Art. No.			
	935115		935116	

One-piece Implant System








SICmax onepiece Screw Implants		
Diameter	 Ø 2.8 mm	
Material: Titanium Grade 4		
Ø Prosthetic connection [mm]	3.0	
Length [mm]	9.5	13.0
Geometry of the connection [mm]	cone with ball Ø 1.8	
	Art. No.	
	935290	935291

SICmax onepiece Accessories



Geometry of the connection: Cone with Torx	 SIC Standard Abutment, straight, for cementation	 SIC Standard Abutment, 15° angle, for cementation	 SIC Wax-Up Abutment, straight, residue-free burn-out	 SIC Retention Attachment, for cementation
Ø Implant connection [mm]	3.0			
Ø Prosthetics [mm]	5.0	4.5	4.6	–
Total height [mm]	6.5	7.5	9.0	3.8
Mod. height [mm]	4.7	5.8	8.0	–
Gingival height [mm]	1.0/1.8		1.0	2.3
Material	Titanium Grade 5		PMMA	Titanium Grade 5 and TiN
	Art. No.			
	935410	935411	935409*	935412**

***Please note: This item can be used as provisional abutment, as well.**

****Compatible with Matrix Set 935717 (see p. 51)**



Art. No.	Art. No.
937111  SIC Insertion Tool, Angle Piece, for SICmax onepiece Length [mm] 20.0 Material stainless steel for surgical devices Geometry of the connection Cone with Torx	935415  SIC Replacement O-Ring for SICmax onepiece, 10 pieces, hard, green Material Silicone, 70 Shore Geometry of the connection [mm] for ball Ø 1.8
935413  SIC Transfer Cap, reposition, for SICmax onepiece Specific height/length [mm] 5.5 Material POM Geometry of the connection Cone with Torx	935416  SIC Replacement O-Ring for SICmax onepiece, 10 pieces, soft, red Material Silicone, 50 Shore Geometry of the connection [mm] for ball Ø 1.8
935408  SIC Lab Implant for SICmax onepiece Diameter prosthetic connection [mm] 3.0 Diameter [mm] 2.8 Total height [mm] 18.0 Material Titanium Grade 4 Geometry of the connection [mm] cone with ball Ø 1.8	936233  SIC Scan Adapter for SICmax onepiece, for digital modelling of CAD/CAM Customized Milled Abutments Precondition: Open CAD/CAM system SimedaCAD, 3 Shape, ExoCAD (Fraunhofer), DentalWings or LaserDenta with output of unencrypted STL-data Manufacturing of Abutments: simeda medical, Luxemburg Geometry of the connection Cone with Torx
935414  SIC O-Ring Attachment for SICmax onepiece Diameter [mm] 4.8 Total height [mm] 3.5 Material Titanium Grade 5 Geometry of the connection [mm] for ball Ø 1.8	

Provisional Implants

SIC Provisional Implants			2 pieces, incl. choice of 2 Abutments, (O-Ring/Conical Abutment)		
Diameter	Ø 2.3 mm		Ø 2.3 mm		
Material: Titanium Grade 4					
Ø Prosthetic connection [mm]	2.75				
Geometry of the connection [mm]	External Square and Ball Ø 1.8				
Length [mm]	10.0		13.0		
	Art. No.				
	935401		935400		

SIC Accessories Provisional Implants

SIC Drills for Provisional Implants


Material: stainless steel for surgical devices		
Diameter [mm]	1.5	1.7
Length [mm]	28.0	34.0
Specific height/length [mm]	13.0	22.0
Geometry of the connection	ISO Adapter Angle Piece	
	Art. No.	
	935229	935407*

*only for bone quality D1

Art. No.	
935404	 SIC Insertion Tool for Provisional Implants Length [mm] 19.0 Material stainless steel for surgical devices Geometry of the connection Inner Square
935402	 SIC Conical Abutment for Provisional Implants Diameter prosthetic connection [mm] 2.75 Diameter [mm] 4.3 Total height [mm] 6.5 Angulation [°] 0 Material Titanium Grade 5 Geometry of the connection Inner Square
935403	 SIC O-Ring Attachment for Provisional Implants Diameter prosthetic connection [mm] 1.8 Diameter [mm] 4.8 Total height [mm] 3.5 Material Titanium Grade 5 Geometry of the connection Ball/O-ring


Prosthetics

SICace



Prosthetics


SICmax



Prosthetics

SIC Gingiva Shapers	p. 28
SIC Impression Technique	p. 30
SIC Fixation Screws	p. 32
SIC Insertion Tools	p. 32
SIC Lab Implants	p. 32
SIC P2F Abutment Dr. Galip Gurel	p. 33
SIC Abutments	p. 34
SIC Multi-Unit Abutment System	p. 41

SIC Bar and Bridge Abutments	p. 45
SIC CAD/CAM	p. 47
Locator®	p. 50
SIC Ball and Socket Attachments	p. 52
Steco Magnetic Attachments	p. 53
Accessories SICmax onepiece	p. 54
Accessories SIC Provisional Implants	p. 55
SIC Prosthetic Tools	p. 56



SIC Surgical / Prosthetic Concept

The inner geometry of the implant is designed as a **precision inner hex** for fitting the abutment components and other system components. The hex features long guide surfaces and the ultimate in manufacturing precision. Another feature is the comparatively long retention screw with a diameter of 1.6 mm. The high, flexible initial tension of the screw reliably prevents it loosening and, combined with the quality characteristics of the hex, ensures reliable retention of the abutment components with outstanding long-term stability. The uniform tightening torque for all retention screws is 20 Ncm.

2-part SIC implants have **platform switching** in the form of a 45° angled, conical implant shoulder. The abutment components are fitted in the prosthetic implant connection diameter.






The assignment between implant and prosthetic diameter is illustrated in the overview on the left.

Material data sheets for all prosthetic abutments are available for downloading in the internet at www.sic-invent.com.

Gingiva Shapers

SIC Gingiva Shapers Ø 3.3 mm – cylindrical

Diameter implant connection 3.3 mm

Compatible with	SICace Ø 3.4 mm SICmax Ø 3.7 mm			SICace Ø 4.0 mm SICmax Ø 4.2 mm	
Material: Titanium Grade 5					
Ø Prosthetics [mm]	3.3				
Gingival height [mm]	2.0	3.0	4.0	5.0	7.0
Art. No.					
	935023	935024	935026	935025	935027*

*Please note: This gingiva shaper should mainly be used for bite registration.

SIC Gingiva Shapers Ø 3.3 mm – concave, anterior

Diameter implant connection 3.3 mm

Compatible with	SICace Ø 3.4 mm SICmax Ø 3.7 mm			SICace Ø 4.0 mm SICmax Ø 4.2 mm	
Material: Titanium Grade 5					
Ø Prosthetics [mm]	4.3				
Gingival height [mm]	1.0	2.0	3.0	4.0	5.0
Art. No.					
	935061	935062	935063	935065	935064






SIC Gingiva Shapers Ø 3.3 mm – concave, posterior

Diameter implant connection 3.3 mm

Compatible with	SICace Ø 3.4 mm SICmax Ø 3.7 mm			SICace Ø 4.0 mm SICmax Ø 4.2 mm	
Material: Titanium Grade 5					
Ø Prosthetics [mm]	5.3				
Gingival height [mm]	1.5	3.0	4.0	5.0	
Art. No.					
	935080	935081	935087	935082	

SIC Gingiva Shapers Ø 4.2 mm – cylindrical

Diameter implant connection 4.2 mm

Compatible with	SICace Ø 4.5 mm SICmax Ø 4.7 mm			SICace Ø 5.0 mm SICmax Ø 5.2 mm	
Material: Titanium Grade 5					
Ø Prosthetics [mm]	4.2				
Gingival height [mm]	2.0	3.0	4.0	5.0	7.0
Art. No.					
	935074	935075	935073	935079	935071*

*Please note: This gingiva shaper should mainly be used for bite registration.

SIC Gingiva Shapers Ø 4.2 mm – concave, anterior

Diameter implant connection 4.2 mm





Compatible with	SICace Ø 4.5 mm SICmax Ø 4.7 mm			SICace Ø 5.0 mm SICmax Ø 5.2 mm	
Material: Titanium Grade 5					
Ø Prosthetics [mm]	5.3				
Gingival height [mm]	1.0	2.0	3.0	4.0	7.0
Art. No.					
	935069	935076	935077	935072	935078
					935070

SIC Gingiva Shapers Ø 4.2 mm – concave, posterior


Diameter implant connection 4.2 mm





Compatible with	SICace Ø 4.5 mm SICmax Ø 4.7 mm			SICace Ø 5.0 mm SICmax Ø 5.2 mm	
Material: Titanium Grade 5					
Ø Prosthetics [mm]	6.3				
Gingival height [mm]	1.5	3.0	4.0	5.0	7.0
Art. No.					
	935088	935083	935089	935084	935090

Impression Technique


SIC Impression Technique		
SIC Transfer Abutment Ø 3.3 mm, reposition (incl. SIC Standard Fixation Screw and Transfer Cap Ø 3.3 mm)		
Compatible with	 SICace Ø 3.4 mm SICmax Ø 3.7 mm	 SICace Ø 4.0 mm SICmax Ø 4.2 mm
Material: Titanium Grade 5		
Ø Implant connection [mm]	3.3	
Gingival height [mm]	4.0	7.0
Specific height/length [mm]	8.8	11.8
Geometry of the connection [mm]	Hex 2.3	
	Art. No.	
	936134	936033

Accessories

SIC Transfer Caps Ø 3.3 mm, 5 pieces Material: POM	
	
	Art. No.
	936030





SIC Impression Technique		
SIC Transfer Abutment Ø 4.2 mm, reposition (incl. SIC Standard Fixation Screw and Transfer Cap Ø 4.2 mm)		
Compatible with	 SICace Ø 4.5 mm SICmax Ø 4.7 mm	 SICace Ø 5.0 mm SICmax Ø 5.2 mm
Material: Titanium Grade 5		
Ø Implant connection [mm]	4.2	
Gingival height [mm]	4.0	7.0
Specific height/length [mm]	8.8	11.8
Geometry of the connection [mm]	Hex 2.3	
	Art. No.	
	936023	936042

Accessories

SIC Transfer Caps Ø 4.2 mm, 5 pieces Material: POM	
	
	Art. No.
	936031





SIC Impression Technique

SIC Transfer Abutment Ø 3.3 mm,
Open Tray Technique (incl. SIC Fixation Screw)








Compatible with	SICace Ø 3.4 mm SICmax Ø 3.7 mm		SICace Ø 4.0 mm SICmax Ø 4.2 mm	
Material: Titanium Grade 5				
	incl. SIC Fixation Screw 17.5 mm	incl. SIC Fixation Screw 22.5 mm	incl. SIC Fixation Screw 22.5 mm	incl. SIC Fixation Screw 27.5 mm
Ø Implant connection [mm]	3.3			
Specific height/length [mm]	9.8		13.8	
Geometry of the connection [mm]	Hex 2.3			
	Art. No.			
	936201	936207	936211	936212

SIC Impression Technique

SIC Transfer Abutment Ø 4.2 mm,
Open Tray Technique (incl. SIC Fixation Screw)


Compatible with	SICace Ø 4.5 mm SICmax Ø 4.7 mm		SICace Ø 5.0 mm SICmax Ø 5.2 mm	
Material: Titanium Grade 5				
	incl. SIC Fixation Screw 17.5 mm	incl. SIC Fixation Screw 22.5 mm	incl. SIC Fixation Screw 22.5 mm	incl. SIC Fixation Screw 27.5 mm
Ø Implant connection [mm]	4.2			
Specific height/length [mm]	9.8		13.8	
Geometry of the connection [mm]	Hex 2.3			
	Art. No.			
	936203	936208	936216	936217

SIC Fixation Screws for Open Tray Technique







Compatible with	 SICace Ø 3.4 mm SICmax Ø 3.7 mm	 SICace Ø 4.0 mm SICmax Ø 4.2 mm	 SICace Ø 4.5 mm SICmax Ø 4.7 mm	 SICace Ø 5.0 mm SICmax Ø 5.2 mm
Material: Titanium Grade 5				
Total height [mm]	17.5	22.5	27.5	
Specific height/length [mm]	3.0	8.0	13.0	
Geometry of the connection [mm]	Torx Ø 2.8 and Inner Hex 1.22			
	Art. No.			
	936525	936527	936526	

SIC Insertion Tool

(manual use)

Compatible	for SIC Fixation Screw, for Transfer Abutment, Open Tray Technique
Material: stainless steel for surgical devices	
Diameter [mm]	7.5
Total height [mm]	6.0
Geometry of the connection [mm]	Torx Ø 2.8
	Art. No.
	937042

SIC Lab Implants

Compatible with	 	 
Material: stainless steel for surgical devices		
Diameter [mm]	3.4	5.0
Ø Prosthetics [mm]	3.3	4.2
Total height [mm]	11.5	
Geometry of the connection [mm]	Hex 2.3	
	Art. No.	
	936133	936232





Abutments

SIC P2F Abutment Dr. Galip Gurel

SIC P2F Abutment Dr. Galip Gurel

Ø 3.3 mm, round (incl. SIC Standard Fixation Screw)







Compatible with	 	
Material: PEEK OPTIMA [®]		
Ø Prosthetics [mm]	3.3	
Gingival height [mm]	2.0	4.0
Art. No.		
	935750* [NEW]	935751* [NEW]

SIC P2F Abutment Dr. Galip Gurel

Ø 3.3 mm, triangle (incl. SIC Standard Fixation Screw)







Compatible with	 	
Material: PEEK OPTIMA [®]		
Ø Prosthetics [mm]	3.3	
Gingival height [mm]	2.0	4.0
Art. No.		
	935752* [NEW]	935753* [NEW]

SIC P2F Abutment Dr. Galip Gurel

Ø 3.3 mm, oval (incl. SIC Standard Fixation Screw)









Compatible with	 	
Material: PEEK OPTIMA [®]		
Ø Prosthetics [mm]	3.3	
Gingival height [mm]	2.0	4.0
Art. No.		
	935754* [NEW]	935755* [NEW]



Abutments, incl. Standard Fixation Screw






SIC Standard Abutments, anterior, straight

Ø 3.3 mm
(incl. SIC Standard Fixation Screw)

Compatible with	<div>  SICace Ø 3.4 mm SICmax Ø 3.7 mm  SICace Ø 4.0 mm SICmax Ø 4.2 mm </div>			
Material: Titanium Grade 5				
Ø Implant connection [mm]	3.3			
Ø Prosthetics [mm]	4.4			
Mod. height [mm]	7.0			
Specific height [mm]	8.0	9.0	10.0	12.0
Gingival height [mm]	1.0	2.0	3.0	5.0
Geometry of the connection [mm]	Hex 2.3			
	Art. No.			
	936163	935800	936153	935801



SIC Standard Abutments, anterior, 15° angle

Ø 3.3 mm
(incl. SIC Standard Fixation Screw)

Compatible with	<div><div> SICace Ø 3.4 mm SICmax Ø 3.7 mm</div><div> SICace Ø 4.0 mm SICmax Ø 4.2 mm</div></div>		
Material: Titanium Grade 5			
Ø Implant connection [mm]	3.3		
Ø Prosthetics [mm]	4.4		
Mod. height [mm]	7.0		
Specific height [mm]	8.0	10.0	12.0
Gingival height [mm]	1.0	3.0	5.0
Geometry of the connection [mm]	Hex 2.3		
	Art. No.		
	936164	936154	936182

SIC Standard Abutments, anterior, 25° angle




Ø 3.3 mm
(incl. SIC Standard Fixation Screw, short)

Compatible with	SICace Ø 3.4 mm SICmax Ø 3.7 mm	SICace Ø 4.0 mm SICmax Ø 4.2 mm
Material: Titanium Grade 5		
Ø Implant connection [mm]	3.3	
Ø Prosthetics [mm]	4.4	
Mod. height [mm]	7.0	
Specific height [mm]	8.0	10.0
Gingival height [mm]	1.0	3.0
Geometry of the connection [mm]	Hex 2.3	
	Art. No.	
	936300 [NEW]	936301 [NEW]

Caution: With SIC Standard Abutment 25° always use SIC Standard Fixation Screw, short




SIC Standard Abutments, posterior, straight

Ø 3.3 mm
(incl. SIC Standard Fixation Screw)

Compatible with	SICace Ø 3.4 mm SICmax Ø 3.7 mm		SICace Ø 4.0 mm SICmax Ø 4.2 mm	
Material: Titanium Grade 5				
Ø Implant connection [mm]	3.3			
Ø Prosthetics [mm]	5.4			
Mod. height [mm]	7.0			
Specific height [mm]	8.5	10.0	12.0	
Gingival height [mm]	1.5	3.0	5.0	
Geometry of the connection [mm]	Hex 2.3			
	Art. No.			
	935727	935728	935802	






SIC Standard Abutments, posterior, 15° angle

Ø 3.3 mm
(incl. SIC Standard Fixation Screw)

Compatible with	SICace Ø 3.4 mm SICmax Ø 3.7 mm			SICace Ø 4.0 mm SICmax Ø 4.2 mm		
Material: Titanium Grade 5						
Ø Implant connection [mm]	3.3					
Ø Prosthetics [mm]	5.4					
Mod. height [mm]	7.0					
Specific height [mm]	8.5		10.0		12.0	
Gingival height [mm]	1.5		3.0		5.0	
Geometry of the connection [mm]	Hex 2.3					
	Art. No.					
	935731		935732		936183	





SIC Standard Abutments, anterior, straight

Ø 4.2 mm
(incl. SIC Standard Fixation Screw)

Compatible with	SICace Ø 4.5 mm SICmax Ø 4.7 mm		SICace Ø 5.0 mm SICmax Ø 5.2 mm		
Material: Titanium Grade 5					
Ø Implant connection [mm]	4.2				
Ø Prosthetics [mm]	5.4				
Mod. height [mm]	7.0				
Specific height [mm]	8.0	9.0	10.0	12.0	14.0
Gingival height [mm]	1.0	2.0	3.0	5.0	7.0
Geometry of the connection [mm]	Hex 2.3				
	Art. No.				
	936165	935803	936152	935804	935805



SIC Standard Abutments, anterior, 15° angle

Ø 4.2 mm
(incl. SIC Standard Fixation Screw)

Compatible with	SICace Ø 4.5 mm SICmax Ø 4.7 mm		SICace Ø 5.0 mm SICmax Ø 5.2 mm	
Material: Titanium Grade 5				
Ø Implant connection [mm]	4.2			
Ø Prosthetics [mm]	5.4			
Mod. height [mm]	7.0			
Specific height [mm]	8.0	10.0	12.0	14.0
Gingival height [mm]	1.0	3.0	5.0	7.0
Geometry of the connection [mm]	Hex 2.3			
	Art. No.			
	936166	936657	936184	936185

SIC Standard Abutments, anterior, 25° angle





Ø 4.2 mm
(incl. SIC Standard Fixation Screw, short)

Compatible with	SICace Ø 4.5 mm SICmax Ø 4.7 mm	SICace Ø 5.0 mm SICmax Ø 5.2 mm
Material: Titanium Grade 5		
Ø Implant connection [mm]	4.2	
Ø Prosthetics [mm]	5.4	
Mod. height [mm]	7.0	
Specific height [mm]	8.0	10.0
Gingival height [mm]	1.0	3.0
Geometry of the connection [mm]	Hex 2.3	
	Art. No.	
	936304 [NEW]	936305 [NEW]

Caution: With SIC Standard Abutment 25° always use SIC Standard Fixation Screw, short





SIC Standard Abutments, posterior, straight

Ø 4.2 mm
(incl. SIC Standard Fixation Screw)

Compatible with	SICace Ø 4.5 mm SICmax Ø 4.7 mm		SICace Ø 5.0 mm SICmax Ø 5.2 mm	
Material: Titanium Grade 5				
Ø Implant connection [mm]	4.2			
Ø Prosthetics [mm]	6.4			
Mod. height [mm]	7.0			
Specific height [mm]	8.5	10.0	12.0	14.0
Gingival height [mm]	1.5	3.0	5.0	7.0
Geometry of the connection [mm]	Hex 2.3			
	Art. No.			
	935729	935730	935726	935806









SIC Standard Abutments, posterior, 15° angle

Ø 4.2 mm
(incl. SIC Standard Fixation Screw)

Compatible with	SICace Ø 4.5 mm SICmax Ø 4.7 mm		SICace Ø 5.0 mm SICmax Ø 5.2 mm	
Material: Titanium Grade 5				
Ø Implant connection [mm]	4.2			
Ø Prosthetics [mm]	6.4			
Mod. height [mm]	7.0			
Specific height [mm]	8.5	10.0	12.0	14.0
Gingival height [mm]	1.5	3.0	5.0	7.0
Geometry of the connection [mm]	Hex 2.3			
	Art. No.			
	935733	935734	936186	936187









SIC "Flex Star" Universal Cast-to Abutment

for non-precious alloy
(incl. SIC Standard Fixation Screw)

Compatible with	 SICace Ø 3.4 mm SICmax Ø 3.7 mm	 SICace Ø 4.0 mm SICmax Ø 4.2 mm	 SICace Ø 4.5 mm SICmax Ø 4.7 mm	 SICace Ø 5.0 mm SICmax Ø 5.2 mm
Material: Ptlr		 for bridgework		 for bridgework
Ø Implant connection [mm]	3.3		4.2	
Ø Prosthetics [mm]	3.6		4.4	
Specific height [mm]	8.8			
Height base [mm]	0.8			
Geometry of the connection [mm]	Hex 2.3	Cone 2.3	Hex 2.3	Cone 2.3
	Art. No.			
	936615	936613	936616	936614









SIC "Flex Star" Universal Cast-to Abutment

high-fusing alloy
(incl. SIC Standard Fixation Screw)








Compatible with	 SICace Ø 3.4 mm SICmax Ø 3.7 mm	 SICace Ø 4.0 mm SICmax Ø 4.2 mm	 SICace Ø 4.5 mm SICmax Ø 4.7 mm	 SICace Ø 5.0 mm SICmax Ø 5.2 mm
Material: AuPt		 for bridgework		 for bridgework
Ø Implant connection [mm]	3.3		4.2	
Ø Prosthetics [mm]	3.6		4.4	
Specific height [mm]	8.8			
Height base [mm]	0.8			
Geometry of the connection [mm]	Hex 2.3	Cone 2.3	Hex 2.3	Cone 2.3
	Art. No.			
	936647	936645	936665	936646

SIC "White Star" Zirconium Oxide Abutment

burn-onable
(incl. SIC Standard Fixation Screw)

Compatible with	 SICace Ø 3.4 mm SICmax Ø 3.7 mm	 SICace Ø 4.0 mm SICmax Ø 4.2 mm	 SICace Ø 4.5 mm SICmax Ø 4.7 mm	 SICace Ø 5.0 mm SICmax Ø 5.2 mm
Material: ZrO ₂ -TZP Colour: white				
Ø Implant connection [mm]	3.3		4.2	
Ø Prosthetics [mm]	3.6		4.8	
Mod. height [mm]	4.0		5.0	
Specific height [mm]	5.5	7.0	6.5	8.0
Gingival height [mm]	1.5	3.0	1.5	3.0
Geometry of the connection [mm]	Hex 2.3			
	Art. No.			
	936148	936149	936150	936151

SIC Accessories for Abutments

Art. No.	   	Art. No.
936658	SIC Standard Fixation Screw Diameter thread [mm] 1.6 Total height [mm] 9.2 Specific height/length [mm] 3.5 Material Titanium Grade 5 Geometry of the connection [mm] 1.6 and Inner Hex 1.22 Compatible with all SIC Standard Abutments	936664  SIC Sleeve for Horizontal Screw (high-fusing alloy) Total height [mm] 1.9 Width spec. [mm] 2.0 Material AuPt Compatible with all SIC Standard Abutments
936659	 SIC Horizontal Screw, length-adjustable (for transversal screw connection) Total height [mm] 5.7 Specific height/length [mm] 2.8 Material Titanium Grade 5 Geometry of the connection [mm] 1.4 and Inner Hex 1.22 Compatible with all SIC Standard Abutments, Titanium	936110  SIC Universal Abutment-Holder Total height [mm] 60.0 Width [mm] 8.0 Material stainless steel for surgical devices Geometry of the connection [mm] Inner Hex 2.3 Compatible with all SIC Standard Abutments


SIC Multi-Unit Abutment System

SIC "Safe on Four"®









The SIC "Safe on Four" system is a further development of the "bar and bridge abutments" system. In the "Safe on Four" system, the bar and bridge abutments and the "Safe on Four" angled standard abutments are directly screw-retained with the respective implant. In this way, a fixed transgingival platform is created over which all further prosthetic and laboratory technical measures are completed. The system is indicated for fixed or removable bridge or full restorations with the stipulation that the distal implants can have a maximum implantation angle of 30°.

The maximum bone availability is utilised distally by displacement of the most distally placed angled implant. During development of the system, particular value was placed on the greatest possible stability of individual components. The straight bar and bridge abutments consist of a two-part design which includes an abutment component with hex and a "Safe on Four" fixation post that, using the long screw shank, ensures maximum continuous loading capacity. The "Safe on Four" universal fixation screws also have a reinforced screw thread.

For further information about the surgical procedure of our Multi-Unit Abutment System please refer to the brochure.














SIC Bar and Bridge Abutment "Safe on Four" (incl. SIC Fixation Post "Safe on Four")






Compatible with	 SICace Ø 3.4 mm SICmax Ø 3.7 mm	 SICace Ø 4.0 mm SICmax Ø 4.2 mm	 SICace Ø 4.5 mm SICmax Ø 4.7 mm	 SICace Ø 5.0 mm SICmax Ø 5.2 mm
Material: Titanium Grade 5				
Ø Implant connection [mm]	3.3		4.2	
Ø Prosthetics [mm]	5.0			
Specific height/length [mm]	4.9	6.4	4.9	6.4
Gingival height [mm]	1.5	3.0	1.5	3.0
Geometry of the connection [mm]	Hex 2.3			
	Art. No.			
	936253	936254	936255	936256

Caution: Tightening torque for SIC Fixation Post "Safe on Four": 30 Ncm with any implant insertion tool SW 2.3









SIC Standard Abutment "Safe on Four" (incl. SIC Fixation Screw, short)

Compatible with	 			 			
Material: Titanium Grade 5							
Ø Implant connection [mm]	3.3			4.2			
Angulation [°]	16	16	30	16	16	30	30
Ø Prosthetics [mm]	5.0						
Specific height/length [mm]	4.9	6.4	4.9	4.9	6.4	4.9	6.4
Gingival height [mm]	1.5	3.0	1.5	1.5	3.0	1.5	3.0
Geometry of the connection [mm]	Hex 2.3						
	Art. No.						
	936259	936261	936260	936263	936262	936264	936265

SIC Accessories for Multi-Unit Abutment System

Art No.		Art No.	
936252	 SIC Gingiva Shaper "Safe on Four", cylindrical Total height [mm] 5.0 Diameter [mm] 5.0 Material Titanium Grade 5 Geometry of the connection Cone "Safe on Four"	936541	 SIC Fixation Screw "Safe on Four" 16.0 mm, for Transfer Abutment, Open Tray Technique Total height (mm) 16.0 Specific height (mm) 6.0 Material Titanium Grade 5 Geometry of the connection (mm) 1.6 and inner hex 1.22
936250	 SIC Transfer Abutment "Safe on Four", reposition (incl. SIC Fixation Screw "Safe on Four" and Transfer Cap Ø 3.3 mm) Total height [mm] 9.5 Diameter [mm] 5.0 Material Titanium Grade 5 Geometry of the connection Cone "Safe on Four"	936251	 SIC Lab Implant "Safe on Four" Total height [mm] 16.0 Diameter [mm] 5.0 Material stainless steel for surgical devices Geometry of the connection Cone "Safe on Four"
936275 [NEW]	 SIC Transfer Abutment "Safe on Four", Open Tray Technique (incl. SIC Fixation Screw 16.0 mm) Total height [mm] 11.0 Diameter [mm] 5.0 Material Titanium Grade 5 Geometry of the connection Cone "Safe on Four"		

SIC Crown Base "Safe on Four" (incl. SIC Fixation Screw "Safe on Four")

Compatible with	 SICace Ø 3.4 mm SICmax Ø 3.7 mm	 SICace Ø 4.0 mm SICmax Ø 4.2 mm	 SICace Ø 4.5 mm SICmax Ø 4.7 mm	 SICace Ø 5.0 mm SICmax Ø 5.2 mm
	 Titanium	 residue-free burn-out	 high-fusing alloy	 for non-precious alloy
Total height [mm]	9.0	10.0	12.4	12.4
Material	Titanium Grade 5	PMMA	AuPt	PtIr
Geometry of the connection	Cone "Safe on Four"			
	Art. No.			
	936270 [NEW]	936276 [NEW]	936272	936273

Art. No.	
936274	SIC Bar Coping "Safe on Four" (incl. SIC Fixation Screw "Safe on Four") Total height [mm] 6.2 Material AuAgCuPt Geometry of the connection Cone "Safe on Four"
936271	SIC Scan Adapter for SIC "Safe on Four" Total height [mm] 6.0 Material PEEK Geometry of the connection Cone "Safe on Four"
936280	SIC Positioning Aid 16° and 30° for "Safe on Four" Material Titanium Grade 5 Geometry of the connection [mm] Hex 2.3
936281	SIC Planning Guide for SIC "Safe on Four" Material Titanium Grade 5









Art. No.	
936257	SIC Fixation Post "Safe on Four", GH 1.5 mm Total height [mm] 9.7 Material Titanium Grade 5 Geometry of the connection Cone "Safe on Four" and Hex 2.3 Compatible with 936253 and 936255
936258	SIC Fixation Post "Safe on Four", GH 3.0 mm Total height [mm] 11.2 Material Titanium Grade 5 Geometry of the connection Cone "Safe on Four" and Hex 2.3 Compatible with 936254 and 936256
936540	SIC Fixation Screw "Safe on Four" Diameter thread [mm] 1.6 Total height [mm] 5.0 Material Titanium Grade 5 Geometry of the connection [mm] 1.6 and Inner Hex 1.22
936529	SIC Fixation Screw, short Total height [mm] 5.0 Material Titanium Grade 5 Geometry of the connection [mm] 1.6 and Inner Hex 1.22
935230	SIC Cutter Ø 5.2 mm for SIC "Safe on Four", with Guidance Tip Diameter [mm] 5.2 Length [mm] 31.0 Material stainless steel for surgical devices Geometry of the connection ISO Adapter Angle Piece





SIC Bar and Bridge Abutment

The conventional “bar and bridge abutment” system features a two-part design consisting of a crown base or bar coping which is screw-retained on the bar and bridge abutment using a central retention screw. The system is suitable for use with fixed or removable bridges or full restorations. The divergence between the implants can be a maximum of 15°.





SIC Bar and Bridge Abutments

Compatible with	 SICace Ø 3.4 mm SICmax Ø 3.7 mm	 SICace Ø 4.0 mm SICmax Ø 4.2 mm	 SICace Ø 4.5 mm SICmax Ø 4.7 mm	 SICace Ø 5.0 mm SICmax Ø 5.2 mm
Material: Titanium Grade 5				
Ø Implant connection [mm]	3.3		4.2	
Ø Prosthetics [mm]	4.6		4.6	
Specific height [mm]	3.5	5.0	3.5	5.0
Cone height [mm]	2.0			
Gingival height [mm]	1.5	3.0	1.5	3.0
Geometry of the connection [mm]	Hex 2.3			
	Art. No.			
	936170	936171	936172	936173





SIC Crown Bases and Bar Coping

Compatible with			
Geometry of the connection : conical, round	 SIC Crown Base, Cast-to, for Bar and Bridge Abutments, Ø 3.3/4.2 mm (for non-precious alloy)	 SIC Crown Base, Cast-to, for Bar and Bridge Abutments Ø 3.3/4.2 mm (high-fusing alloy)	 SIC Bar Coping for Bar and Bridge Abutments Ø 3,3/4,2 mm
Ø Prosthetics [mm]	4.7		4.6
Total height [mm]	10.4	10.4	4.5
Material	PtIr	AuPt	AuAgCuPt
Art. No.			
	936617	936670	936680




SIC Fixation Screws

Compatible with	Bar and Bridge Abutments			
Material: Titanium Grade 5 Geometry of the connection [mm]: 1.6 and Inner Hex 1.22	 for Crown Base	 for Crown Base	 for Bar Coping	 for Bar Coping
Length [mm]	12.7	14.2	11.2	12.7
Gingival height [mm]	1.5	3.0	1.5	3.0
Art. No.				
	936430	936431	936432	936433

SIC Soldering Implants for Bar and Bridge Abutment

Compatible with		
Material: stainless steel for surgical devices		
Ø Prosthetics [mm]	4.6	
Length [mm]	19.0	
Geometry of the connection [mm]	conical, round	
	Art. No.	
	936691	936692









SIC Bars

Compatible with		
Material: AuPtPdAgCu	 SIC Standard Bar, round, Ø 1.9 mm, gold	 SIC Bar Clip, gold
Length [mm]	50.0	
	Art. No.	
	936682	936683

CAD/CAM

SIC CAD/CAM







SIC Bonding Base CAD/CAM, straight,
(incl. SIC Standard Fixation Screw)

Compatible with	 SICace Ø 3.4 mm SICmax Ø 3.7 mm	 SICace Ø 4.0 mm SICmax Ø 4.2 mm	 SICace Ø 4.5 mm SICmax Ø 4.7 mm	 SICace Ø 5.0 mm SICmax Ø 5.2 mm
Material: Titanium Grade 5		for bridgework 		for bridgework 
Ø Implant connection [mm]	3.3		4.2	
Ø Prosthetics [mm]	3.5		4.4	
Mod. height [mm]	4.0			
Specific height [mm]	4.3			
Gingival height [mm]	0.3			
Geometry of the connection [mm]	Hex 2.3	Cone 2.3	Hex 2.3	Cone 2.3
	Art. No.			
	936190	936191	936196	936197

Please note: This item can be used as provisional abutment, as well.

SIC CAD/CAM







SIC Bonding Base CAD/CAM, 15° angle,
(incl. SIC Standard Fixation Screw)

Compatible with	 SICace Ø 3.4 mm SICmax Ø 3.7 mm	 SICace Ø 4.0 mm SICmax Ø 4.2 mm	 SICace Ø 4.5 mm SICmax Ø 4.7 mm	 SICace Ø 5.0 mm SICmax Ø 5.2 mm
Material: Titanium Grade 5				
Ø Implant connection [mm]	3.3		4.2	
Ø Prosthetics [mm]	3.5		4.4	
Mod. height [mm]	4.0			
Specific height [mm]	4.3			
Gingival height [mm]	0.3			
Geometry of the connection [mm]	Hex 2.3			
	Art. No.			
	936192		936198	

Please note: This item can be used as provisional abutment, as well.

SIC CAD/CAM







SIC Bonding Base CAD/CAM Ø 3.3 mm, straight, CEREC (incl. SIC Standard Fixation Screw)

Compatible with	 SICace Ø 3.4 mm SICmax Ø 3.7 mm	 SICace Ø 4.0 mm SICmax Ø 4.2 mm	 SICace Ø 4.5 mm SICmax Ø 4.7 mm	 SICace Ø 5.0 mm SICmax Ø 5.2 mm
Material: Titanium Grade 5				
Ø Implant connection [mm]	3.3		4.2	
Ø Prosthetics [mm]	4.1		4.5	
Mod. height [mm]	4.7			
Specific height [mm]	5.0			
Gingival height [mm]	0.3			
Geometry of the connection [mm]	Hex 2.3			
	Art. No.			
	936188		936189	

Please note: This item can be used as provisional abutment, as well.

SIC CAD/CAM







SIC Wax-Up Base CAD/CAM Ø 3.3 mm (incl. SIC Standard Fixation Screw)

Compatible with	 SICace Ø 3.4 mm SICmax Ø 3.7 mm	 SICace Ø 4.0 mm SICmax Ø 4.2 mm	 SICace Ø 4.5 mm SICmax Ø 4.7 mm	 SICace Ø 5.0 mm SICmax Ø 5.2 mm
Material: Titanium Grade 5				
Ø Implant connection [mm]	3.3		4.2	
Ø Prosthetics [mm]	4.3			
Mod. height [mm]	7.5			
Specific height [mm]	8.3			
Gingival height [mm]	0.8			
Geometry of the connection	Hex 2.3			
	Art. No.			
	936701		936703	

Please note: This item can be used as provisional abutment, as well.



SIC Scan Adapter

for digital modelling of CAD/CAM Customized Milled Abutments
(incl. SIC Standard Fixation Screw)







Compatible with	 SICace Ø 3.4 mm SICmax Ø 3.7 mm	 SICace Ø 4.0 mm SICmax Ø 4.2 mm	 SICace Ø 4.5 mm SICmax Ø 4.7 mm	 SICace Ø 5.0 mm SICmax Ø 5.2 mm
				
Ø Implant connection [mm]	3.3		4.2	
Geometry of the connection [mm]	Hex 2.3			
	Art. No.			
	936234		936235	

Precondition: Open CAD/CAM system
SimedaCAD, 3 Shape, ExoCAD (Fraunhofer),
DentalWings or LaserDenta with output of
unencrypted STL-data
Manufacturing of Abutments:
simeda medical, Luxemburg







SIC Accessories for CAD/CAM

Art. No.	
936528	SIC Standard Fixation Screw red, for Ceramic Abutment CAD/CAM  <div> Length [mm] 9.1 Head diameter [mm] 2.2 Material Titanium Grade 5 Geometry of the connection [mm] 1.6 and Inner Hex 1.22 Compatible with SIC CAD/CAM ceramic abutments </div>
936529	SIC Fixation Screw, short  <div> Length [mm] 7.3 Head diameter [mm] 1.9 Material Titanium Grade 5 Geometry of the connection [mm] 1.6 and Inner Hex 1.22 Compatible with SIC Bonding Base angled and Standard Abutment "Safe on Four" and 25° Abutments </div>

Locator®

Locator® Attachments						
incl. Locator® Matrix Set 935717						
Compatible with	SICace Ø 3.4 mm SICmax Ø 3.7 mm			SICace Ø 4.0 mm SICmax Ø 4.2 mm		
Material: Titanium Grade 5 and TIN						
Ø Implant connection [mm]	3.3					
Gingival height [mm]	0.0	1.0	2.0	3.0	4.0	5.0
Art. No.						
	935702	935710	935703	935711	935704*	935708*

*Indicated with a minimum implant diameter of 4,0 mm or 4 locator abutments on one framework






Locator® Attachments						
incl. Locator® Matrix Set 935717						
Compatible with	SICace Ø 4.5 mm SICmax Ø 4.7 mm			SICace Ø 5.0 mm SICmax Ø 5.2 mm		
Material: Titanium Grade 5 and TIN						
Ø Implant connection [mm]	4.2					
Gingival height [mm]	0.0	1.0	2.0	3.0	4.0	5.0
Art. No.						
	935705	935712	935706	935713	935707	935709











Locator[®] Replacement Male 4 pieces each

Compatible with	all Locator [®] Attachments						
							
Material	Nylon						
Retention force [kg]	2.3	1.4	0.7	no retention	1.8	0.9	0.5
Colour	clear	pink	blue	grey	green	orange	red
	Art. No.						
	935718	935719	935724	935723	935720	935725	935721






Locator[®] Accessories

Art. Nr.					Art. Nr.				
935717		Locator[®] Matrix Set, 5-part			935722		Locator[®] Insertion Tool, Angle Piece		
		Material	Nylon				Material	stainless steel for surgical devices	
		Compatible with	all Locator [®] Attachments				Compatible with	all Locator [®] Attachments	
935715		Locator[®] Impression Post			935714		Locator[®] Tool, 3-part		
		Material	Titanium Grade 4				Material	stainless steel for surgical devices	
		Compatible with	all Locator [®] Attachments				Compatible with	all Locator [®] Attachments	
935716		Locator[®] Lab Implant							
		Material	Titanium Grade 4						
		Compatible with	all Locator [®] Attachments						

SIC Ball and Socket Attachments

Compatible with	 SICace Ø 3.4 mm SICmax Ø 3.7 mm	 SICace Ø 4.0 mm SICmax Ø 4.2 mm	 SICace Ø 4.5 mm SICmax Ø 4.7 mm	 SICace Ø 5.0 mm SICmax Ø 5.2 mm
Material: Titanium Grade 5				
Ø Implant connection [mm]	3.3		4.2	
Ø Prosthetics [mm]	3.3		4.2	
Gingival height [mm]	2.0	4.0	2.0	4.0
Art. No.				
	936157	936158	936016	936017

Inner Matrices

Compatible with	SIC Ball and Socket Attachments				
					
	Inner Assembly Matrix, blue (for model)	Inner Matrix, gold, suitable for activation	Inner Matrix, yellow, light	Inner Matrix, green, standard	Inner Matrix, red, strong
Material	Titanium Grade 5	AuPt	Titanium Grade 5		
Retention force [kg]	–	0.8	0.6	0.8	1.2
Art. No.					
	936001	936002	936005	936004	936003





Accessories Ball and Socket Attachments

Art. Nr.		Art. No.	
936008	Retention Cap, universal Material Titanium Grade 5 Compatible with SIC Ball and Socket Attachments	936015	SIC Paralleling Tool Material stainless steel for surgical devices Compatible with SIC Ball and Socket Attachments
936013	Spacer for Retention Cap, universal Material POM Compatible with SIC Ball and Socket Attachments	937040	Activator/Deactivator for Inner Matrix "gold" Material stainless steel for surgical devices Compatible with SIC Ball and Socket Attachments
936610	Lab Implant for Ball and Socket Attachments Material stainless steel for surgical devices Compatible with SIC Ball and Socket Attachments	936006	Torque Ratchet, 4 Ncm, Matrix Material stainless steel for surgical devices Compatible with SIC Ball and Socket Attachments
936010	Polymerisation Aid (PCV ring), 10 pieces Material PVC Compatible with SIC Ball and Socket Attachments		
936014	Block-Out Aid (distance plate), 2 pcs. Material Zn Compatible with SIC Ball and Socket Attachments		

Steco Magnetic Attachment








Art. No.		Art. No.	
935095	Steco Magnetic Abutment Ø 4.2 mm Material Titanium Grade 4	935098	Steco Lab Implant for Magnetic Abutment Material Titanium Grade 4
935096	Steco Magnet for Prosthesis Material Titanium Grade 4	935099	Steco Insertion Tool for Magnetic Abutment Material Titanium Grade 4
935097	Steco Positioning Cuff for Magnetic Abutment Material Dental Silicone		

SICmax onepiece Accessories

Geometry of the connection: Cone with Torx	 SIC Standard Abutment, straight, for cementation	 SIC Standard Abutment, 15° angle, for cementation	 SIC Wax-Up Abutment, straight, residue-free burn-out	 SIC Retention Attachment, for cementation
Ø Implant connection [mm]	3.0			
Ø Prosthetics [mm]	5.0	4.5	4.6	–
Total height [mm]	6.5	7.5	9.0	3.8
Mod. height [mm]	4.7	5.8	8.0	–
Gingival height [mm]	1.0/1.8		1.0	2.3
Material	Titanium Grade 5		PMMA	Titanium Grade 5 and TiN
	Art. No.			
	935410	935411	935409*	935412**



*Please note: This item can be used as provisional abutment, as well.

**Compatible with Matrix Set 935717 (see p. 51)

Art. No.		Art. No.	
	937111 SIC Insertion Tool, Angle Piece, for SICmax onepiece Length [mm] 20.0 Material stainless steel for surgical devices Geometry of the connection Cone with Torx		935415 SIC Replacement O-Ring for SICmax onepiece, 10 pieces, hard, green Material Silicone, 70 Shore Geometry of the connection [mm] for ball Ø 1.8
	935413 SIC Transfer Cap, reposition, for SICmax onepiece Specific height/length [mm] 5.5 Material POM Geometry of the connection Cone with Torx		935416 SIC Replacement O-Ring for SICmax onepiece, 10 pieces, soft, red Material Silicone, 50 Shore Geometry of the connection [mm] for ball Ø 1.8
	935408 SIC Lab Implant for SICmax onepiece Diameter prosthetic connection [mm] 3.0 Diameter [mm] 2.8 Total height [mm] 18.0 Material Titanium Grade 4 Geometry of the connection [mm] cone with ball Ø 1.8		936233 SIC Scan Adapter for SICmax onepiece, for digital modelling of CAD/CAM Customized Milled Abutments Precondition: Open CAD/CAM system SimedaCAD, 3 Shape, ExoCAD (Fraunhofer), DeantalWings or LaserDenta with output of unencrypted STL-data Manufacturing of Abutments: simeda medical, Luxemburg Geometry of the connection Cone with Torx
	935414 SIC O-Ring Attachment for SICmax onepiece Diameter [mm] 4.8 Total height [mm] 3.5 Material Titanium Grade 5 Geometry of the connection [mm] for ball Ø 1.8		

SIC Accessories Provisional Implants

SIC Drills for Provisional Implants

Material: stainless steel for surgical devices		
Diameter [mm]	1.5	1.7
Length [mm]	28.0	34.0
Specific height/length [mm]	13.0	22.0
Geometry of the connection	ISO Adapter Angle Piece	
	Art. No.	
	935229	935407*

*only for bone quality D1

Art. No.

935404



SIC Insertion Tool for Provisional Implants

Length [mm] 19.0
Material stainless steel for surgical devices
Geometry of the connection Inner Square

935402



SIC Conical Abutment for Provisional Implants

Diameter prosthetic connection [mm] 2.75
Diameter [mm] 4.3
Total height [mm] 6.5
Angulation [°] 0
Material Titanium Grade 5
Geometry of the connection Inner Square

935403





SIC O-Ring Attachment for Provisional Implants

Diameter prosthetic connection [mm] 1.8
Diameter [mm] 4.8
Total height [mm] 3.5
Material Titanium Grade 5
Geometry of the connection Ball/O-ring


Prosthetic Tools

SIC Prosthetic Tray


		
	SIC Prosthetic Tray TR, equipped (incl. TR Screwdrivers – short and long – and TR Torque Ratchet incl. Adapter)	SIC Prosthetic Tray, empty
Length/Height/Width [mm]	155/25/55	
	Art. No.	
	935526	935523

Art. No.	<div><div></div><div></div><div></div><div></div></div>	
937127	SIC Torque Ratchet, Titanium, incl. Adapter for Angle Piece Instruments	
	Diameter [mm]	7.0
	Length [mm]	99.0
	Width [mm]	12.0
	Material	Titanium Grade 5
937108	SIC TR Adapter for Angle Piece Instruments	
	Length [mm]	19.5
	Material	stainless steel for surgical devices
	Compatible with	SIC Angle Piece Instruments
936110	SIC Universal Abutment-Holder	
	Total height [mm]	60.0
	Width [mm]	8.0
	Material	stainless steel for surgical devices
	Geometry of the connection [mm]	Inner Hex 2.3
	Compatible with	all SIC Standard Abutments




SIC Screwdriver, Angle Piece

Material: stainless steel for surgical devices		
	short, hex. 1.2 mm	long, hex. 1.2 mm
Specific height/length [mm]	5.0	15.0
	Art. No.	
	937032	937031

SIC Insertion Tool (manual use)

Compatible	for Fixation Screw, for Transfer Abutment, Open Tray Technique
Material: stainless steel for surgical devices	
Diameter [mm]	7.5
Total height [mm]	6.0
Geometry of the connection [mm]	Torx Ø 2.8
	Art. No.
	937042

SIC TR Screwdriver

Material: stainless steel for surgical devices			
	extra short, hex. 1.2 mm	short, hex. 1.2 mm	long, hex. 1.2 mm
Length [mm]	17.0	25.0	32.5
Geometry of the connection [mm]	Hex 1.2		
	Art. No.		
	937130	937128	937129

Standard Surgery

SIC Bone Condenser

The SIC Bone Condenser, with an instrument design based on an idea suggested by Dr. A. Weidmann, enables atraumatic condensing of the implant site in the cancellous bone. Bone preparation with the new, patented instrument geometry of the condenser attachment greatly improves primary stability during implant placement in soft bone.



SIC Titanium Ratchet

The SIC titanium ratchet (TR) combines maximum precision, secure handling, improved durability with an attractive modern design. The one-piece ratchet body is made of a titanium alloy and the snap-on ratchet head is made of stainless steel, guaranteeing high protection against corrosion and easy and thorough cleaning, care and maintenance. For checking torque, the ergonomically designed handle has an individually calibrated and scaled cam follower which is employed at torques of up to 45 Ncm. When used as a surgical ratchet, without torque control, operation without the cam follower allows transmission of a maximum torque of 80 Ncm. The TR adapters engage in the ratchet seating and lock into it securely. The adapter for angled instruments (article no. 937108) has a special "Hex socket" for transmitting torques of up to 80 Ncm.



SIC Surgical Tray	p. 60
SIC Drill Systems	p. 60
Marking Drill	p. 60
Short Drills	p. 61
Drills for Depth Stop	p. 61
SIC Depth Stops for Drills	p. 62
SIC Brackets for Depth Stops	p. 63
SIC Crestal Drills	p. 63
SIC Reamers	p. 64
SIC Bone Taps	p. 64

SIC Bone Condenser	p. 65
SIC Surgical Accessories	p. 66
SIC Insertion Tool Surgery	p. 68

SIC Drilling System with a Depth Stop

The SIC drilling system with depth stop is a very flexible system for reliable, quick implant placement. The drills can be used with or without the depth stop. A secure friction fit guarantees the high functionality of the depth stop. This is achieved by a clip on the drill depth stop which can simply be attached over the working section of the drill.





SIC Surgical Tray



The SIC Surgical Tray is characterised by its optimal arrangement and ergonomic design. Due to its small dimensions, it fits into the smallest sterilising machines even when it is completely filled. The number of instruments is reduced to the necessities. The drill system can be plugged into the tray in modular fashion.

There is the option of adding drill depth stops to the tray.



All three SIC invent implant lines are compatible with the instruments of the SIC surgical tray. With this "SIC one4all Concept™", we meet the demands of the dental practitioner for ease of use and cost-effectiveness.



Surgical Trays

SIC Surgical Tray		
Compatible with	SICace SICmax SICmax onepiece	
		
	SIC Surgical Tray with Titanium Ratchet, equipped	SIC Surgical Tray universal, empty
Length/Height/Width [mm]	250/25/150	
	Art. No.	
	933105	935510

Drill Systems

SIC Marking-Drill	
Compatible with all SIC implant diameters	
Material: stainless steel for surgical devices	
Diameter [mm]	1.7
Length [mm]	34.0
Specific height/length [mm]	14.5
Geometry of the connection	ISO Adapter Angle Piece
	Art. No.
	935194

SIC Short Drills

Compatible with						
Material: stainless steel for surgical devices		 "Smart Drill"				
Diameter [mm]	2.0	2.0/2.5/2.8	3.1	3.25	3.75	4.25
Length [mm]	25.5					
Specific height/length [mm]	11.5					
Geometry of the connection	ISO Adapter Angle Piece					
	Art. No.					
	935222	935223	935224	935225	935226	935228








Please note: Only to be used with implant lengths up to 11.5 mm.

SIC Drills for Depth Stop

Compatible with						
Material: stainless steel for surgical devices		 "Smart Drill"				
Diameter [mm]	2.0	2.0/2.5/2.8	3.1	3.25	3.75	4.25
Length [mm]	34.0					
Specific height/length [mm]	14.5					
Geometry of the connection	ISO Adapter Angle Piece					
	Art. No.					
	935214	935215	935216	935217	935220	935218

SIC Depth Stops for Drills








up to Ø 3.25 mm

To use with implant diameter						
Material: stainless steel for surgical devices						
Diameter [mm]	4.7					
for L ¹ * [mm]	6.0	7.5	9.5	11.5	13.0	14.5
Specific height/length [mm]	13.9	12.4	10.4	8.4	6.9	5.4
Compatible with	Depth Stop Drills up to Ø 3.25 mm					
Geometry of the connection	ISO Adapter Angle Piece					
	Art. No.					
	935245	935246	935247	935248	935249	935250

*L¹ = Implant Length



SIC Depth Stops for Drills

up to Ø 4.25 mm



To use with implant diameter						
Material: stainless steel for surgical devices						
Diameter [mm]	5.7					
for L ¹ [mm]	6.0	7.5	9.5	11.5	13.0	14.5
Specific height/length [mm]	13.9	12.4	10.4	8.4	6.9	5.4
Compatible with	Depth Stop Drills up to Ø 4.25 mm					
Geometry of the connection	ISO Adapter Angle Piece					
	Art. No.					
	935256	935251	935252	935253	935254	935255

*L¹ = Implant Length









SIC Brackets for Depth Stops

 <p>equipped</p>	 <p>empty</p>
Art. No.	
935512	935511

Compatibility of Depth Stops

Diameter [mm]	2.0	2.0/2.5/2.8	3.1	3.25	3.75	4.25
 Depth stop without marking	935214	935215	935216	935217		
 Depth stop with marking					935220	935218







SIC Crestal Drills

Compatible with	 SICace Ø 3.4 mm SICmax Ø 3.7 mm	 SICace Ø 4.0 mm SICmax Ø 4.2 mm	 SICace Ø 4.5 mm SICmax Ø 4.7 mm	 SICace Ø 5.0 mm SICmax Ø 5.2 mm
Material: stainless steel for surgical devices				
Diameter [mm]	3.3	3.75	4.25	4.75
Length [mm]	31.0			
Geometry of the connection	ISO Adapter Angle Piece			
	Art. No.			
	935187	935192	935167	935193









SIC invent | Standard Surgery

SIC Reamer









for SICace Length 6.0 mm

Compatible with	 Ø 4.0 mm	 Ø 4.5 mm	 Ø 5.0 mm
Material: stainless steel for surgical devices			
Diameter [mm]	3.9	4.4	4.9
Length [mm]	31.0 mm		
Geometry of the connection	ISO Adapter Angle Piece		
	Art. No.		
	935127 [NEW]	935128 [NEW]	935129 [NEW]

SIC Bone Taps

Compatible with	 SICace Ø 3.4 mm SICmax Ø 3.7 mm	 SICace Ø 4.0 mm SICmax Ø 4.2 mm	 SICace Ø 4.5 mm SICmax Ø 4.7 mm	 SICace Ø 5.0 mm SICmax Ø 5.2 mm
Material: stainless steel for surgical devices				
Diameter [mm]	3.4	4.0	4.5	5.0
Length [mm]	34.0			
Specific height/length [mm]	14.5			
Geometry of the connection	ISO Adapter Angle Piece			
	Art. No.			
	935185	935186	935168	935190

SIC Bone Condenser by Dr. A. Weidmann

Compatible with					
Material: stainless steel for surgical devices					
Diameter [mm]	–	2.6	3.2	3.4	4.3
Length [mm]	140.0	27.0			
Specific height/length [mm]	14.5				
	Art. No.				
	937200	937202	937203	937204	937205

SIC Bone Condenser Set by Dr. A. Weidmann




Consisting of Tray, Handle, 4 inserts

Art. No.	
937220	

SIC Bone Condenser





The SIC Bone Condenser, with an instrument design based on an idea suggested by Dr. A. Weidmann, enables atraumatic condensing of the implant site in the cancellous bone. Bone preparation with the new, patented instrument geometry of the condenser attachment greatly improves primary stability during implant placement in soft bone.







The bone condenser is used in ascending sequence until the required final diameter is reached, whereby the sizes, depth marks and colour coding correspond to those of the SIC implant system. By tapping in the bone condenser the local bone is displaced apically and the bone condensed minimally invasively by **rotating the condenser radially (clockwise)**, slowly and gently by at least 90 degrees.






Surgical Accessories

SIC Gingiva Punch "Safety Punch"

Compatible with		
Material: stainless steel for surgical devices		
Diameter [mm]	3.5	4.4
Length [mm]	20.5	
	Art. No.	
	937151	937150

Art. No.	
937323 	SIC Probe for Gingiva Punch, universal Diameter [mm] 1.9 Length [mm] 35.0 Material stainless steel for surgical devices Geometry of the connection [mm] Hex 1.22
Art. No.	
935212 	SIC Drill extension Length [mm] 27.0 Material stainless steel for surgical devices Geometry of the connection ISO Adapter Angle Piece Compatible with SIC Angle Piece Instruments
Art. No.	
935230 	SIC Cutter Ø 5.2 mm for SIC "Safe on Four", with Guidance Tip Diameter [mm] 5.2 Length [mm] 31.0 Material stainless steel for surgical devices Geometry of the connection ISO Adapter Angle Piece

SIC Depth Gauge

Compatible with		
Material: Titanium Grade 5		
Diameter [mm]	2.0	2.8
Length [mm]	23.5	
Specific height/length [mm]	20.0	
	Art. No.	
	935165	935166







SIC Counter Ratchet with "Easy Handle"

Diameter [mm] 1.6
Length [mm] 28.0
Geometry of the connection Hex 3.0
Material Titanium Grade 5



	Art. No.
	937122 [NEW]

SIC Distance Control "Surgery"

Compatible with	 	 
Material: Titanium Grade 5		
Ø Implant connection [mm]	3.3	4.2
Length [mm]	15.5	
Specific height [mm]	10.5	
Geometry of the connection [mm]	cylindrical 2.0	
	Art. No.	
	935320	935322

Art. No.



935351



SIC X-Ray Balls Ø 5.0 mm, 10 pieces




Diameter [mm] 5.0
Material stainless steel for surgical devices

Z9913





SIC Dental Implant Pass, 20 pieces

SIC X-Ray Positioning Devices



			
	SICace®, 1.25:1 and 1:1	SICmax®, 1.25:1 and 1:1	SICmax onepiece, 1.25:1 and 1:1
Art. No.			
	935007	935008	935009

Insertion Tools Surgery




SIC Insertion Tools, Angle Piece

Material: stainless steel for surgical devices		
		long
Length [mm]	20.0	28.0
Geometry of the connection [mm]	Hex 3.0	
	Art. No.	
	937113	937112



SIC TR Insertion Tools






Material: stainless steel for surgical devices		
	"direct", short	"direct", long
Length [mm]	22.5	32.5
Geometry of the connection [mm]	Hex 2.3	
	Art. No.	
	937102	937103

SIC TR Screwdrivers

Material: stainless steel for surgical devices			
	extra short, hex. 1.2 mm	short, hex. 1.2 mm	long, hex. 1.2 mm
Length [mm]	17.0	25.0	32.5
Geometry of the connection [mm]	Hex 1.2		
	Art. No.		
	937130	937128	937129

SIC Screwdriver, Angle Piece

Material: stainless steel for surgical devices		
	short, hex. 1.2 mm	long, hex. 1.2 mm
Specific height/length [mm]	5.0	15.0
	Art. No.	
	937032	937031

Art. No.			Art. No.		
	937127	SIC Torque Ratchet, Titanium, incl. Adapter for Angle Piece Instruments Diameter [mm] 7.0 Length [mm] 99.0 Width [mm] 12.0 Material Titanium Grade 5		937122 [NEW]	SIC Counter Ratchet with "Easy Handle" Diameter [mm] 1.6 Length [mm] 28.0 Geometry of the connection Hex 3.0 Material Titanium Grade 5
	937108	SIC TR Adapter for Angle Piece Instruments Length [mm] 19.5 Material stainless steel for surgical devices Compatible with SIC Angle Piece Instruments		935300	SIC "Easy Screw" designed by Dr. G. Bayer Length [mm] 175.0 Material stainless steel for surgical devices Compatible with ISO Adapter Angle Piece and Implant Insertion Post
	937109	SIC TR Insertion Tool S, short Length [mm] 19.0 Material stainless steel for surgical devices Geometry of the connection [mm] Hex 3.0			

Guided Surgery



SIC Surgical Tray, Guided Surgery	p. 72
SIC GS Gingiva Punches	p. 72
SIC Drill Systems	p. 73
SIC GS Countersinks	p. 73
SIC GS Pilot Drills	p. 73
SIC GS Extension Drills	p. 74
SIC GS Bone Taps	p. 75
SIC GS Insertion Tools	p. 75
SIC GS Drill Keys	p. 76
SIC GS Sleeves	p. 76

■ SIC Guided Surgery

3D diagnosis in combination with prosthetically oriented backwards planning increases reliability for determining the optimal tooth position.



SIC Guided Surgery is a software-independent surgical system for template-guided, navigated implant insertion. Important characteristics are compactness, efficiency and ergonomics of the instrumentarium. Maximum flexibility due to open connection to current planning tools, variability due to the possibility of laboratory or industrial production of the guiding templates, surgical freedom with maximum functionality and precision predominated during the conception and development of the system.

- Software-independent, open instrument set
- Implemented in the planning tools: SimPlant[®] (Materialise Dental), CeHa imPLANT[®] (med 3D), coDiagnostiX[®] (Straumann[®]), SKYplanX (bredent), SICAT Implant (SICAT GmbH & Co. KG), smop Planning Solution (Swissmeda AG), Nemotec (Software Nemotec, S.L.), Implant Studio (3Shape).
- Fabrication of the guide centrally by Materialise Dental, SICAT GmbH & Co. KG or in a local dental laboratory
- Guidance of implant placement using the guide template
- Maximum flexibility for the operator (no fixed depth stops)
- Master sleeve Ø 5.2 mm for standard indications and Master sleeve Ø 3.1 mm for lateral and lower incisors





Guided Surgery









SIC Surgical Tray, Guided Surgery

Compatible with	<div><div></div>SICace Ø 3.4 mm SICmax Ø 3.7 mm</div> <div><div></div>SICace Ø 4.0 mm SICmax Ø 4.2 mm</div> <div><div></div>SICace Ø 4.5 mm SICmax Ø 4.7 mm</div> <div><div></div>SICace Ø 5.0 mm SICmax Ø 5.2 mm</div>
	<div> equipped</div> <div> empty</div>
Length/Height/Width [mm]	250/25/150
	Art. No.
	933110935518






SIC GS Gingiva Punches

Compatible with	<div><div></div><div></div></div>	<div><div></div><div></div></div>
Material: stainless steel for surgical devices		
Diameter [mm]	4.2	5.2
Length [mm]	26.0	
Specific height/length [mm]	10.0	
Geometry of the connection	ISO Adapter Angle Piece	
	Art. No.	
	937154	937155






SIC GS Countersinks

Compatible with	 SICace Ø 3.4 mm SICmax Ø 3.7 mm	 SICace Ø 4.0 mm SICmax Ø 4.2 mm	 SICace Ø 4.5 mm SICmax Ø 4.7 mm	 SICace Ø 5.0 mm SICmax Ø 5.2 mm
Material: stainless steel for surgical devices				
Diameter [mm]	3.4	3.9	4.4	4.9
Length [mm]	29.5			
Specific height/length [mm]	15.5			
Geometry of the connection	ISO Adapter Angle Piece			
	Art. No.			
	935550	935551	935566	935552










SIC GS Pilot Drill

Compatible with	   
Material: stainless steel for surgical devices	
Diameter [mm]	2.0
Length [mm]	35.5
Specific height/length [mm]	21.5
Geometry of the connection	ISO Adapter Angle Piece
	Art. No.
	935553










SIC GS Pilot Drill, long

Compatible with	   
Material: stainless steel for surgical devices	
Diameter [mm]	2.0
Length [mm]	40.0
Specific height/length [mm]	26.0
Geometry of the connection	ISO Adapter Angle Piece
	Art. No.
	935554









SIC GS Extension Drills

Compatible with					
Material: stainless steel for surgical devices					
	"Smart Drill"				
Diameter [mm]	2.0/2.5/2.8	3.1	3.25	3.75	4.25
Length [mm]	35.5				
Specific height/length [mm]	21.5				
Geometry of the connection	ISO Adapter Angle Piece				
	Art. No.				
	935555	935556	935557	935567	935558







SIC GS Extension Drills, long

Compatible with					
Material: stainless steel for surgical devices					
	"Smart Drill"				
Diameter [mm]	2.0/2.5/2.8	3.1	3.25	3.75	4.25
Length [mm]	40.0				
Specific height/length [mm]	26.0				
Geometry of the connection	ISO Adapter Angle Piece				
	Art. No.				
	935559	935560	935561	935568	935562



SIC GS Bone Taps

Compatible with	 SICace Ø 3.4 mm SICmax Ø 3.7 mm	 SICace Ø 4.0 mm SICmax Ø 4.2 mm	 SICace Ø 4.5 mm SICmax Ø 4.7 mm	 SICace Ø 5.0 mm SICmax Ø 5.2 mm
Material: stainless steel for surgical devices				
Diameter [mm]	3.4	4.0	4.5	5.0
Length [mm]	35.5			
Geometry of the connection	ISO Adapter Angle Piece			
	Art. No.			
	935563	935564	935569	935565





SIC GS Insertion Tools, Angle Piece


Compatible with	   	
Material: stainless steel for surgical devices		 long
Length [mm]	23.0	28.0
Geometry of the connection [mm]	Hex 3.0	
	Art. No.	
	937115	937114

SIC GS Drill Keys

Compatible with	(see Surgical Guideline SIC Guided Surgery)			
Material: stainless steel for surgical devices	 and for Bone Tap for Sleeve Ø 5.2 mm			
Diameter 1 [mm]	2.0, Sleeve Ø 3.1 mm	2.0	3.10	3.75
Diameter 2 [mm]	2.4 (Tap), Sleeve Ø 5.2 mm	2.8	3.25	4.25
For sleeve Ø [mm]	5.2			
	Art. No.			
	935580	935581	935582	935583

SIC GS Sleeves

Compatible with: (see Surgical Guideline SIC Guided Surgery)		
Material: Titanium Grade 4		
Diameter [mm]	3.1	5.2
Outer diameter [mm]	4.0	6.0
Art. No.		
	935590	935591




Planning dimensions for
sleeve positioning in the
guide template



Art. No.	
935592	SIC GS Fixation Post, for Sleeve Ø 5.2
	Diameter [mm] 5.2 Specific height [mm] 9.0

SIC Guided Surgery – Overview of Tools

Standard Surgery



Guided Surgery



Guided Surgery

Biomaterial

SIC nature graft

SIC nature graft is a purely biological, anorganic augmentation material of phycogene origin.

SIC nature graft is distinguished by an interconnecting porosity, a honey-comb-like structure and a coarse surface. These characteristics guarantee optimum osteoconduction and rapid bone regeneration.



Bone Grafting

p. 78

SIC nature graft

p. 79

Bone Grafting Material

SIC nature graft



SIC nature graft 1×0.5 ml
(granule size 0.3–1.0 mm)





SIC nature graft 1×1.0 ml
(granule size 0.3–1.0 mm)




SIC nature graft 1×2.0 ml
(granule size 0.3–1.0 mm)

Material	hydroxyapatite of phycogenic origin		
	Art. No.		
	510808	510816	510824


SIC Certificates



CERTIFICATE



This is to certify that the company



SIC invent AG
Birmannsgasse 3
4055 Basel
Schweiz

with the organizational units/sites as listed in the annex


has implemented and maintains a **Quality Management System**.

Scope:
Development, manufacturing and distribution of dental implant systems, risk class I instruments
and bioresorbable bone filling material

Through an audit, documented in a report, performed by DQS Medizinprodukte GmbH, it was
verified that the management system fulfills the requirements of the following standard:


EN ISO 13485 : 2012 + AC : 2012

Certificate registration No.	293238 MP2012
Certificate unique ID	170592665
Effective date	2014-04-17
Expiry date	2016-12-31
Frankfurt am Main	2014-04-17




Deutsche
Akkreditierungsstelle
D-ZM-16021-01-01

DQS Medizinprodukte GmbH



Frank Gralchen
Managing Director



Dr. Thomas Feldmann
Head of Certification Body

August-Schanz-Straße 21, 60433 Frankfurt am Main, Tel. +49 (0) 69 95427-263, medical.devices@dqs.de

1 / 2

510(k) approved by FDA, USA



EC-CERTIFICATE

(Full quality assurance system)



This is to certify that the company



SIC invent AG

Birmannsgasse 3
4055 Basel
Schweiz

has implemented and maintains a full quality assurance system which applies to the products at every stage from design to final controls.

Through an audit, documented in a report, performed by DQS Medizinprodukte GmbH, it was verified that the management system fulfills the requirements of

Annex II – excluding Section 4 of Council Directive 93/42/EEC concerning medical devices

with respect to the following medical devices:

Dental Implants, Abuments, Accessories as listed in the annex

The manufacturer is subject to surveillance according to Annex II, Section 5. The CE marking with the Notified Body Identification Number (0297) may be affixed on the devices listed in the certificate. An EC Design Examination Certificate according to Annex II, Section 4 is required for class III devices covered by this certificate. The certificate is in the case of class I(s) devices (I(s) = class I products placed on the market in sterile conditions) limited to the aspects of manufacture concerned with securing and maintaining sterile conditions. The certificate is in the case of class I(m) devices (I(m) = class I devices with a measuring function) limited to the aspects of manufacture concerned with the conformity of the products with the metrological requirements.

Certificate registration No.	293238 MR2
Certificate unique ID	170592329
Effective date	2014-04-08
Expiry date	2018-12-31
Frankfurt am Main	2014-04-08

DQS Medizinprodukte GmbH

Frank Graichen
Managing Director

Dr. Thomas Feldmann
Head of Certification Body

August-Schanz-Straße 21, 60433 Frankfurt am Main, Tel. +49 (0) 69 95427-253, medical.devices@dqs.de

DQS Medizinprodukte GmbH is a Notified Body according to Council Directive 93/42/EEC concerning medical devices with the Identification Number 0297.

1 / 2

20130805/02/EN/30308



SIC invent AG

Birmannsgasse 3
CH-4055 Basel
Tel.: +41 (0)61 260 24 60
Fax: +41 (0)61 261 39 68
contact.switzerland@sic-invent.com

SIC invent Deutschland GmbH

Bühlstraße 21
D-37073 Göttingen
Tel.: 0800 742 46 83 68
Fax: 0800 74 23 29 38
contact.germany@sic-invent.com

SIC invent Austria GmbH

Kohlmarkt 7/Stg. 2/58
A-1010 Wien
Tel.: +43 (0)1 533 70 60
Fax: +43 (0)1 533 70 60 50
contact.austria@sic-invent.com

SIC invent Asia-Pacific Inc.

6F, Hyojin Bldg, 1540-5 Seocho-dong
Seocho-gu, Seoul, 137-070, Korea (ROK)
Tel.: +82 2 585 9700
Fax: +82 2 584 4411
contact.korea@sic-invent.com

SIC invent North America Inc.

contact.usa@sic-invent.com

www.sic-invent.com
www.my-sic.com