A.G.M Dental Material

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Polycarboxylate Cement

Introduction:

AGM Zinc polycarboxylate cement is a dental cement with excellent adhesion which is used for cementation of veneers, orthodontic appliances, inlays, intraoral splints, partial dentures, permanent metal restorations, temporary restorations as well as coating/base of Cavities in amalgam and composite materials. This cement contains two parts including powder and liquid, the powder part consists of zinc oxide and magnesium oxide and the liquid part contains a polymer aqueous solution based on polyacrylic acid.

Special Features:

AGM polycarboxylate cement offers a strong force /strength to traction and compression. It has excellent adhesion to the surface of enamel and dentin. Its liquid component is polyacrylic acid, which owing to its high molecular weight prevents pulp sensibility since their molecules are large and cannot traverse the peritubular space acting as sealants.

Composition

Powder:

Zinc oxide
Magnesium oxide

Liquid:

Polyacrylic acid

Packaging

25 gm Powder

12 ml Liquid

Instructions for use:

1- Mixing ratio:

For cementation: 1 measure of powder to 3 drop of liquid according to desired consistency.

For linings and basing: 2 measure of powder to 3 drop of liquid.

2- An appropriate ratio of powder and liquid should be placed on a clean and dry mixing pad/glass slab.

3- With suitable spatula incorporate all of the powder into liquid rapidly.

Notes:

-Mixing process should be executed at room temperature.

-The cement should be used after the completion of mix -while the appearance is glossy.

- Avoid adding extra powder and liquid after the mixing has been started

-Setting time is proportional to temperature and humidity.

Shelf Life: 2 years

Physical Feature:

Type Feature	Luting	Base/Line
Compressive strength (1 day)	58 Mp	65 Mp
Film thickness	25 µm	-
Mixing time	40 second	30 second
Working time	90 second	70 second
Setting time	7.5 min	6 min

Use & Storage

Store in a cool and dark place

Warning

- Kept out of reach of children

- Because of adhesive nature of polycarboxylate cement all the instruments comprising spatulas, glass slabs and etc. should be washed with plenty amount of water instantly after usage.

Emergency Measures

Avoid contact with eyes, skin and oral tissues:

- In the case of eyes contact, rinse instantly with plenty amount of water and seek medical attention.

- In the case of skin contact wash with soap and water.

- In the condition of ingestion dilute by drinking water.

Standards:

ISO	13485	9917	29022	10993-5
	10993-3		10993-10	

Andishe Gostar Masoud Co.

Tehran, Iran

AGMDENTALMATERIAL.COM

FDA Number: 45362906

IMD Number: 14400281

Product Code: AGDC05