



Preparation protocol.













Preparation protocol with CITO mini® preparation instruments.

Taking into account different bone qualities.

Adapt preparation protocol according to the indication and situation of each patient.

The insertion of CITO mini® implants in soft bone is not indicated.



		Soft bone quality			Medium bone quality			Hard bone quality			
											
		ø 1.8 ³	ø 2.2	ø 2.5	ø 1.8 ³	ø 2.2	ø 2.5	ø 1.8 ³	ø 2.2	ø 2.5	
Marking drill ²					X	X	X	X	X	X	
Stepped drill for ø 1.8 mm ¹		The insertion of CITO mini® implants in soft bone is not indicated.			min. 6.0 mm			X			
Stepped drill for ø 2.2 mm ¹					min. 6.0 mm				X		
Stepped drill for ø 2.5 mm ¹								min. 6.0 mm			X

¹ The insertion depth/length of the drill depends on the desired primary stability and the bone quality. The depth markings/scaling help to avoid damaging neighboring structures (e.g. nerves) due to perforations that are too deep.

² Optional use.

³ Not intended for maxillary prostheses.