

Essential endo motor

User-friendly

- Torque controlled motor for safe operation
- Broad torque range for operation of several NiTi systems
- Individual selection of torque and speed values for all NiTi systems
- Additional individual program for storage of max. 15 personal settings
- Convenient battery operation, motor can also be used while charging



→ Focus on treatment

- Pre-programmed torque and speed settings for Mtwo®, FlexMaster® and Gates
- Intuitive menu navigation
- Simple operation and clear display



→ Automatic safety functions

- Automatic reverse rotation when set torque limit is reached
- Acoustic warning signals
 - at 75 % of the set torque values
 - when rotating in reverse direction

→ High quality contra-angle

- Small head for good visibility and easy handling
- Fast file change



Range of products

REF 1153 VDW.SILVER®

Including Sirona 6:1 contra-angle and foot control

REF 1079 Separate Sirona 6:1 contra-angle with small head

REF 1155 VDW.SILVER® (REF 1153) + Mtwo® StarterKit

Contents Mtwo® StarterKit:
2 Mtwo® of sizes 10/04, 15/05, 20/06, 25/06 each
1 Mtwo® of sizes 30/05, 35/04, 40/04, 25/07 each
1 Mtwo® SystemBox
3 plastic training blocks
2 x 3ml FileCare® EDTA
1 user CD-ROM and instruction brochure

Åhgrens Dental AB
Vegavägen 28
184 50 Åkersberga
www.ahgrens-dental.se
info@ahgrens-dental.se
Tel: 08-540 884 00

REF 1154 VDW.SILVER® (REF 1153) + FlexMaster® StarterKit

Contents FlexMaster® StarterKit:
Each 1 FlexMaster® BasicKit 21 and 25 mm
1 FlexMaster® SystemBox
3 plastic training blocks
2 x 3ml FileCare® EDTA
1 instruction brochure



VDW GmbH
P.O. Box 830954 • 81709 Munich • Germany
Tel. +49 89 62734-0 • Fax +49 89 62734-304
www.vdw-dental.com • info@vdw-dental.com

Endo Easy Efficient®