



Pole Ideal Pars Co.  
www.medpip.com

## Drying rack

Drying and holding laboratory glassware after washing

P.I.P. drying rack is used to safely dry and store frequently used or fragile labware. it is made of resistant coated metal with 65 pegs.

- Useable for all kinds of beakers, Erlenmeyer flasks, laboratory flasks, tubes and any other laboratory glassware.

P.I.P. drying racks can be used on laboratory sink and mounted on the wall.

GS1 Code	Material	Model	Outer dimensions
6260807503447	Metal	65-pegs	420x240x620



1<sup>st</sup> floor, No.12, Naghdi St., Jahantab St., Motahari Ave, Tehran, Iran

FAX:+98 21 88765561/+98 21 88767159 Tel: 88765847 - 88765915 - 88545922-9

Postal code:1576635714 P.O. Box:15875-9483

@poleideal

@poleideal

+98 912 334 0197

Pole Ideal Pars Co.