

RING COLOR	CLINIC	TIPS	RECOMMENDED POWER
GREEN	Periodontics 	P2L	3
		P2R	3
		TK1-1S	2
		TK1-1L	2
		TK2-1L	2
		TK2-1R	2
		PH1	2
		PH2L	2
		PH2R	2
		H1	2
		H2L	2
		H2R	2
		H3	2
		H4L	2
		H4R	2
YELLOW	Endodontics 	ETBD	7
		ET18D	10
		ET20	7
		ET25	7
		ET25S	7
		ET25L	7
		ET20D	7
		ET40	7
		ET40D	7
		SO4	7
		IRR20/21	9
		IRR20/25	9
		IRR25/21	9
		IRR25/25	9
		K 15/21	6
		K 25/25	6
		CAP1	10
		CAP2	10
		CAP3	10
	Surgical endodontics 	AS3D	7
		AS6D	7
		AS9D	6
		ASLD	7
		ASRD	7
		P14D	7
		P15LD	7
		P15RD	7
	Conservatrice and restorative dentistry 	PM2	7
		PM3	10
		PMS2	7
		PMS3	10
		PMV4	10
		PMV5	10
		PMV6	10
BLUE	Prophylaxis 	1	14
		1S	14
		2	14
		3	14
		10P	14
		10X	12
		10Z	12
	Conservatrice and restorative dentistry 	EX1	12
		EX2	12
		EX3	12
		EXL	12
		EXR	12
		PM1	15
		PMS1	15
		PM4	14
		PMV1	15
		PMV2	15
		PMV3	15
ORANGE	Conservatrice and restorative dentistry 	5AE	20
		ETPR	20